GLOBUS

ITALIAN EXCELLENCE

-  **ELETTROSTIMOLATORI** - Manuale di utilizzo
-  **ELECTROSTIMULATORS** - User Manual
-  **ELECTROESTIMULADORES** - Manual de usuario
-  **ÉLECTROSTIMULATEURS** - Manuel de l'utilisateur

ELITE SII

GENESY SII

DUO TENS

GLOBUS

ITALIAN EXCELLENCE

**MANUALE DI UTILIZZO****pg. 4****USER MANUAL****pg. 45****MANUAL DE USUARIO****pg. 84****MANUEL DE L'UTILISATEUR****pg. 125**

EGREGIO CLIENTE

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DEAR CUSTOMER

**THANK YOU FOR CHOOSING A GLOBUS PRODUCT.
WE REMAIN AT YOUR DISPOSAL FOR ANY ASSISTANCE OR
ADVICE YOU MAY NEED**

DOMINO s.r.l.

via Vittorio Veneto 52

31013 - Codognè - TV - Italy

Tel. (+39) 0438.7933

Fax. (+39) 0438.793363

e-mail: info@globuscorporation.com

www.globuscorporation.com

Table of contents

TECHNICAL FEATURES	47
INTENDED USE	48
EQUIPMENT	49
LABELLING AND SYMBOLS	51
PANEL/KEYBOARD	55
PREPARATION TO THE USE OF THE DEVICE	56
Warnings before use	56
Warnings during the use	57
SIDE EFFECTS AND CONTRAINDICATIONS	58
Contraindications	58
USER GUIDE OF THE DEVICE	59
Programming	60
PROGRAM LIST	62
GENERAL NOTES ON ELECTRODE POSITIONING	67
TREATMENT DESCRIPTIONS	70
SPORT	70
BEAUTY	71
TENS	72
ACTION PRINCIPLES	73
Muscular electrostimulation	73
Electrostimulation advantages	73
Mechanism of the muscular contraction	74
Tens	76
Incontinence programs	77
Warnings	77
Maintenance	77
MAINTENANCE AND CLEANING	78
WARRANTY	80
FREQUENTLY ASKED QUESTIONS	81
EMC accompanying documents	82

TECHNICAL FEATURES

Device

Size: 13 X 8 X 2 cm

Weight: 220 g.

Case: in ABS

Protection level: IP 22

Storage and transportation temperature: from -10°C to 45°C

Max. relative humidity: 30% - 75%

Made in China

Charger

Brand:  Switching Power Adaptor

Model: DK7-065-0200-EU

Input: 100-240V ~ - 50-60Hz - 0,06A

Output: 5,8 V = - 200 mA

Polarity: 

Battery Pack: Ni-MH AAA 4,8V 800 mAh

Use conditions

Temperature: from 0°C to 35°C

Max. relative humidity: from 15% to 93%

Atmospheric pressure: from 700 hPa to 1060 hPa

The values indicate the limits allowed if the product or its accessories are not in the original package.

Technical features of the currents EMS and TENS:

Channels available: Channel 1-2

Constant current: Yes

Intensity: 0-100 mA p.p. per channel

Wave form: Rectangular, biphasic, symmetric, compensated

Working frequency: 1-150 Hz

Recovery frequency: 1-150 Hz

Pulse amplitude: 50-400 μ s

Working time: from 1 to 30 seconds

Recovery time: from 0 to 1 minute

Frequency modulation range: continuous variation from 1 to 150 Hz

Min. modulation time: 3 seconds

Amplitude modulation range: continuous variation from 50 to 400 μ s

Ionophoresis

Channels available: Channels 1

Constant current: Yes

Min. Intensity: 0 mA/1000 Ohm

Max. Intensity: 10 mA/1000 Ohm step. 0.1 mA/1000 Ohm

Min. Time: 1 minute

Max. Time: 99 minutes

Disposal of the device

Do not throw the device or parts of it into the fire; dispose of the product in the specialized centers and respecting the regulations in force in your Country. When the product has to be disposed of, the user can give it back to the retailer when purchasing a new device.

A correct separate waste collection or the compliance with the above-mentioned prescriptions contribute to avoiding possible negative effects on the environment and the health and promote the reuse and/or recycle of the materials of which the device is composed. The illegal disposal of the product entails the application of administrative fines according to applicable regulations.

Declaration of conformity

The device has been manufactured in compliance with applicable technical standards and has been certified, in compliance with Directive 93/42/EEC as amended by directive 2007/47 on medical devices, by the Notified Body Kiwa Cermet Italia, Via Cadriano 23, 40057 Granarolo Dell'Emilia (BO) Italy (n. 0476), in order to ensure product safety.

INTENDED USE

The estimated usable life of the product is 5 years. It is advisable to return the device to the producer and/or authorized center to perform security and maintenance checks every 2 years. The number of treatments depends on the battery charge. The useful life of the battery is estimated at 6 months, after which it is recommended to replace it.

The risk class of the devices:

Elite SII **IIa**

Genesy SII **IIb**


Duo TENS **IIa**

The electrostimulators are designed to be used in the following operating environments:

- domestic environment;
- rehabilitation in general;
- pain treatments in general;
- for aesthetic and sport purposes (CE0476 does not refer to non-medical treatments).

The use of this equipment is allowed to the patient (appropriately informed of the conditions of use of the apparatus).

EQUIPMENT

The electrostimulator is supplied complete of cables and electrodes to use: therefore, opening the package, it is necessary to check that the basic equipment is complete. If some elements should be missing, contact immediately the authorized retailer where you purchased the product.

Control carefully the integrity of the device and its electrodes.

A. Device

B. 2 colored electrode connection cables (for EMS, TENS and IONO treatments)

C. 1 bag containing 4 reusable self-adhesive electrodes (50 x 50 mm)

(It is suggested to use these electrodes for small areas such as upper limbs, calves, cervical areas...)

D. 1 bag containing 4 reusable self-adhesive electrodes (50 x 90 mm)

(It is suggested to use these electrodes for big areas such as thighs, abdomen and glute...)

E. Charger (see technical features)

F. Carrying bag

G. User manual and warranty

Accessories that are not included in the equipment (to be purchased separately)

The device can be used with the optional accessories listed below (it is possible to see their features in the website www.globuscorporation.com). In order to purchase these accessories, please contact your dealer.

REF	Name	Description
G1188	Vaginal probe	Probe for the treatment of incontinence and for the reinforcement and/or relaxation of the pelvic floor.
G0757	Anal probe	Probe for the treatment of incontinence and for the reinforcement and/or relaxation of the pelvic floor.
G1156	Motor point pen	It helps finding the best positioning of the electrodes
G0479	Kit conductive elastic bands for thighs	Kit conductive elastic bands for thighs Bands can be used instead of electrodes, and are recommended for aesthetic and beauty treatments.
G0480	Kit conductive elastic bands for thighs and arms Fitness Top	Kit conductive elastic bands for thighs and arms Bands can be used instead of electrodes, and are recommended for aesthetic and beauty treatments.
G0487	Fast Band	Abdominal band for treatment on abdomen, gluteus and back - 98 cm
G0489	Fast Pad	Special reusable electrodes, specifically suitable for aesthetic treatments on thighs and gluteus.
G0885	Big ionophoresis electrode	Carbon electrode + pouch 60x85 mm
G0439	Kit 2 splitting cables	This accessory is used to split the cables in order to use more electrodes at the same time.
G0451	Face electrodes	Facial electrodes 32 mm with cable
G0869	Globus gel 250 ml	Conductive gel tube of 250 ml

CONNECTIONS

Cable connections

It is possible to use one or two cables according as you have chosen to use one or two channels of the unit.

To connect the cable connector to the device plug it into the appropriate inlet in the upper part of the unit.

The inlets are situated exactly under the correspondent channel (left or right). Insert the cable at the channel you want to use.

For IONOPHORESIS use only the Channel 1.

Electrode application

Take the electrodes from the original package; all new electrodes have a seal on the package. Be sure that the device is off. To start, connect the two cable plugs to the electrodes, then disconnect the electrodes from their position and apply them on the skin. To place the electrodes correctly, see the pictures included in this manual.

After the use, place the electrodes in their original position again.

ATTENTION: do not unplug the electrodes if the unit is working.

LABELLING AND SYMBOLS

	It refers to the manufacturer
	Warning
	This symbol indicates that the device complies with the directives on medical devices (93/42/EEC 47/2007/EEC). The number of the notified body is 0476.
	It indicates that this is a class II device.
	It indicates that this device has BF-type applied parts.
	WEEE symbol (Waste of Electrical and Electronic Equipment). Recycling symbol. The WEEE symbol used for this product indicates that it cannot be treated as a household waste. The proper disposal of the product will contribute to protecting the environment. For further information on the recycling of this product, please contact the concerned office of your local body, the household waste management company or the store where the product was purchased.
	It indicates that the product has been designed in compliance with the directive 2011/65/EEC

	It informs the operator that s/he must read the manual before using the device.
	It informs the client of the compulsory conduct
	It indicates the ideal temperature for the storage and transportation of the product.
IP22	Informes the resistance to water ingress
SN	It indicates the serial number of the device.
	It refers to the expiry date of the product
LOT	It refers to the production lot
	It refers to the manufacturing date.
	It indicates the pressure of the environment in which the device and the accessories are transported and stored.
	It indicates the humidity of the environment where the device and its accessories are used and stored
	Indoor use

Devices

Model:  **ELITE SII**  

Electrostimulator Type: A3R
Max Power: 4.8Vdc 150mA
Output: 100mA (1000ohm)
Battery Pack: 800mA 4.8Vdc

Power Input: [®]
Model: DK7-065-0200-EU
Input: 100-240Vac 50-60Hz 0.06A
Output: 5.8Vdc 200mA

SN *A3RWWYYXXXXX*

   0476

 **DOMINO s.r.l.** - Via Vittorio Veneto, 52 31013 Codognè (TV) - IT
Elettrostimolatore - Electrostimulator - Electroestimulador
Electrostimulateur - Elektostimulator

Model:  **GENESY SII**  


Electrostimulator Type: A3R
Max Power: 4.8Vdc 150mA
Output: 100mA (1000ohm)
Battery Pack: 800mA 4.8Vdc

Power Input: [®]
Model: DK7-065-0200-EU
Input: 100-240Vac 50-60Hz 0.06A
Output: 5.8Vdc 200mA

SN *A3RWWYYXXXXX*


   0476

 **DOMINO s.r.l.** - Via Vittorio Veneto, 52 31013 Codognè (TV) - IT
Elettrostimolatore - Electrostimulator - Electroestimulador
Electrostimulateur - Elektostimulator

**SN**


The first 4 digits of the serial number indicate the week and year of manufacture of the device you purchased (for example, if the code is *** 2319 *****, it means that the device was manufactured in week 23 of 2019).

Charger


Electrodes

	It indicates the dimensions of the electrode
	It indicates the number of electrodes contained in the package
REF	It indicates the product code
CE	It refers to product certification and indicates that it complies with directive 2001/95/EC updated as 2014/357/EU
	It indicates the storage temperature of the electrodes


- Clean and degrease the skin.
- Do not apply the electrode on wounds or injured skin.
- Connect the cable connector to the electrode connector.
- Remove the electrode.
- Apply on the skin.
- Start the program.
- At the end, turn off and put the electrode back in the package.
- Electrodes are for personal use.
- Do not remove the electrode by grabbing the connector.
- Electrodes should not touch each other.
- Do not apply the electrodes on the temples, the neck and in a transthoracic way.
- Do not leave the electrodes in the car.

Display and interface

PANEL/KEYBOARD

ON/OFF/OK BUTTON	It confirms the selection. It pauses the running program. 3" = On/Off.
LEFT/BACK Button	To scroll the selection left. To go back to the previous selection. 3" = It returns to the previous phase while the program is running.
P+ Button	It moves the selection upwards. It increases the intensity of the 4 channels simultaneously, while a program is running
P- Button	It moves the selection downwards. It increases the intensity of the 4 channels simultaneously, while a program is running
RIGHT/USER button:	It moves the selection to the right. 3" = It moves to the next phase, while a program is running
* Button	It starts and stops the contraction while "Action Now" programs are running (where "Action Now" programs are included).
Fn (Runtime) button	If pressed together with other buttons, it permits the user to modify their function. Moreover, it selects Runtime mode (editing time, frequency and amplitude)
Intensity button	It increases/decreases the stimulation intensity of the corresponding channel.

PREPARATION TO THE USE OF THE DEVICE

Mandatory behavior

For maximum safety, the device must be used following the rules and the limitations of the user manual.

The manufacturer declines all responsibility with reference to a different use from what is indicated in this manual.

The full or partial reproduction in any form and by any electronic or mechanical means of the texts and/or pictures contained in this manual without the written authorization of the manufacturer is forbidden.

Treatments should not be performed on skin lesions.

If the package, the cable or the connector of the power supply show signs of wear or damage, replace it instantly.

The device should be connected to the mains with its power supply; before doing so, make sure that the power system complies with the directives in force in your country. Make sure that the power supply will be easily unplugged.

Warnings before use

We suggest avoiding the use of the device together with other electronic devices, especially with those used to maintain vital functions. For a correct use of the Electromedical device, make reference to the tables attached. If it is necessary to use the device near or together with other devices, pay attention to its working.

- It is recommended to read carefully the entire operating manual before using the device; keep carefully this operating manual;
- The device can emit current value above 10mA_{rms};
- Before each use always check the integrity of the unit. This is a fundamental requirement for carrying out therapies, do not use the unit if either the buttons or cables are defective or malfunctioning;
- It must be used only by people over 18 who are able to understand and take action; - It should not be used for purposes other than neuromuscular stimulation as described in this manual;
- It must be used following the indications and under the physician or physiotherapist's control;
- It must be used with the supplied electrodes intended for transcutaneous neuromuscular stimulation;
- It must be kept out of the reach of children;
- ECG monitoring devices may not operate properly when electrostimulation is working; - Do not use the device in trans-thoracic modality because it may cause cardiac arrhythmia, imposing its frequency over the heart one. Do not stimulate the pectoral and dorsal muscles simultaneously;

- If there is any health problem, it must be used only after seeing a doctor;
- A simultaneous connection of the patient to a high frequency electrosurgery device can cause burns where the electrodes are placed and the electrostimulator can result damaged;
- After turning on the device, make sure that in the display the software version and the model of the device appear, it means that the device is working and it is ready to be used;

If not, or in the display all the segments do not appear, turn it off and on again. If the problem persists, contact the service center and do not use the device.

An unexpected switching-off of the device means that the battery has run-down. Charge it according to the instructions in the paragraph: how to charge the batteries.

Warnings during the use

While using the electrostimulator some warnings should be followed:

- In case of damaged cables, they must be replaced with original parts and no longer used;
- Use only the electrodes supplied by the company;
- The device must be kept out of the reach of any pet that could damage it and contaminate it with parasites;
- Pay particular attention when the current density for every electrode is above 2mA/cm² (effective value);
- The cables of the electrostimulation should not be wrapped around people's neck to avoid any risk of strangulation and suffocation;
- The mobile and fixed radio communications devices could influence the functioning of the electromedical device: see the tables attached to this manual.

Preventative measures to take using the device for incontinence treatments.

- Patients with extra-urethral incontinence should not be treated with the electrostimulator;
- Patients suffering excessive incontinence due to evacuation problems should not be treated with the electrostimulator;
- Patients with severe urinary retention to the upper urinary passages should not be treated with the stimulator;
- Patients with total peripheral denervation of the pelvic floor should not be treated with the stimulator;
- Patients suffering of a total/subtotal prolapse of the uterus/vagina should be stimulated with extreme care;
- Patients with infections to the urinary passages should be treated for these symptoms before starting the stimulation treatment;
- Before removing or touching the probe it is necessary to turn off the stimulator or to regulate the intensity of both the channels on 0,0 mA;

- The treatment is a personalized medical prescription: do not borrow the stimulator to other persons.

SIDE EFFECTS AND CONTRAINDICATIONS

Isolated incidents of skin irritation may occur in subjects with high epidermal sensitiveness. In case of an allergic reaction to the electrode gel, suspend the treatment and contact a specialist.

If during the treatment signs of tachycardia and extrasystole appear, suspend the treatment and contact your physician.

Contraindications

Do not use the device in the following cases:

- Stimulation of the front part of the neck (carotid sinus);
- Pacemaker weavers;
- Patients with tumor diseases (see your oncologist);
- Stimulation of the brain region;
- Pains whose etiology is unknown;
- Sores and dermatological diseases;
- Severe traumas;
- Stimulation on recent scars;
- Pregnancy;
- It is strictly forbidden to use the electrostimulator in the ocular area;
- Near body areas with metallic implants or intratissue metals (prostheses, osteosynthetic devices, coils, screws, orthopedic plates), when using monophasic current, interferential, or continuous current, ionophoresis.

It is recommended to use the device with caution in people with capillary fragility, as an excessive stimulation could cause capillary ruptures.


For those programs that consist in **anal or vaginal stimulation**, in addition to the general warnings, the following contraindications have to be considered:

- Presence of sexually transmitted diseases.
- Presence of total prolapsed uterus.
- Presence of urinary infections.
- Presence of severe or chronic dermatological diseases.
- Presence of ischemic tissues, wounds, skin or damaged or irritated mucosa, and/or in presence of infections.
- Presence of contraceptive vaginal ring.
- After an invasive or ablative surgical intervention, which has not been totally cured.
- Patients subject to bleedings or using anticoagulants.


- In case of a weak immune system due to an immunosuppressive disease or while using immunosuppressive medicines.
- Presence of overflow incontinence
- Presence of obliterated urethra.
- In case of atrophy of the involved orifice that may cause tissue injuries.

USER GUIDE OF THE DEVICE

Start

Program selection and modification

Programming

By pressing the Fn button for 2 seconds in the starting screen, you will have the possibility to program the device through 6 ways (PED1 – PED6). To move from a parameter to the next one press OK. For all programs it is possible to change the total duration from 5 ' to 100' in steps of 5 '. Default duration 20 '.

PED1:

It is possible to set the values of frequency and duration of the pulse.

From 1 to 150 Hz in steps of 1 Hz.

From 50 to 250 μs in steps of 10 μs .

Default value: 100 Hz / 200 μs .

PED2:

It is possible to regulate the pulse duration with the frequency set on 100 Hz. From 50 to 250 μs in steps of 10 μs .

Default value: 200 μs .

PED 3:

It is possible to set the frequency values.

From 1 to 150 Hz in steps of 1 Hz.

Default values: 100 Hz

It is possible to modulate the pulse duration from 100 μs to 250 μs .

PED 4:

It is possible to set the values of frequency and pulse duration.


From 1 to 150 Hz in steps of 1

Hz.

From 50 to 320 μs in steps of 10 μs .

The following values are not editable: $t_{\text{UP}} = 1 \text{ s}$; $t_{\text{HOLD}} = 5 \text{ s}$; $t_{\text{DOWN}} = 1 \text{ s}$; $t_{\text{OFF}} = 5 \text{ s}$

Default values: 30 Hz / 250 μs .

**PED 5:**

It is possible to set the frequency values.

From 1 to 150 Hz in steps of 1 Hz.

Default values: 30 Hz

It is possible to modulate the pulse duration from 100 μs to 320 μs .

PED 6:

It is possible to set the values of t_{HOLD} e t_{OFF} .

From 1 to 30s in steps of 1s.


It is possible to set the frequency values.

From 1 to 150 Hz in steps of 1 Hz.

Value of the pulse duration set on 300 μ s.

Default value $t_{HOLD} = 2$ s; $t_{OFF} = 10$ s; 30 Hz.

The values of t_{UP} and t_{DOWN} are not modifiable and have a duration of 1s.

PROGRAM LIST

Name	Type	Duo Tens	Elite S2	Genesy S2
DEMO	SPORT		P 1	P 1
Motor point pen	SPORT		P 2	P 2
Warm-up lower limbs	SPORT		P 3	P 3
Pre-competition warm-up lower limbs	SPORT		P 4	P 4
Maximum strength lower limbs	SPORT		P 5	P 5
Endurance strength lower limbs	SPORT		P 6	P 6
Explosive strength lower limbs	SPORT		P 7	P 7
Reactivity lower limbs	SPORT		P 8	P 8
Capillarization lower limbs	SPORT		P 9	P 9
Active recovery lower limbs	SPORT		P 10	P 10
Muscle cool down lower limbs	SPORT		P 11	P11
Endurance lower limbs	SPORT		P 12	
Decontracting lower limbs	SPORT		P 13	P12
Decontracting	SPORT	P 1		
Abdominal maximum strength	SPORT		P 14	

Maximum strength	SPORT	P 2		
Abdominal endurance strength	SPORT		P 15	
Abdominal active recovery	SPORT		P 16	
Abdominal cool down	SPORT		P 17	
Cool down	SPORT	P 3		
Abdominal capillaries	SPORT		P 18	
Abdominal decontracting	SPORT		P 19	
Warm-up upper limbs	SPORT		P 20	
Pre-competition warm-up upper limbs	SPORT		P 21	
Maximum strength upper limbs	SPORT		P 22	
Endurance strength upper limbs	SPORT		P 23	
Explosive strength upper limbs	SPORT		P 24	
Reactivity upper limbs	SPORT		P 25	
Capillarization upper limbs	SPORT		P 26	
Active recovery upper limbs	SPORT		P 27	
Muscle cool down upper limbs	SPORT		P 28	
Endurance upper limbs	SPORT		P 29	
Decontracting upper limbs	SPORT		P 30	
Leg anti-cellulite lipolysis	BEAUTY		P 31	P13
Hip-gluteus anti-cellulite lipolysis	BEAUTY		P 32	
Abdominal area anti-cellulite lipolysis	BEAUTY		P 33	
Upper limb anti-cellulite lipolysis	BEAUTY		P 34	
Leg drainage	BEAUTY		P 35	P14
Hip-gluteus drainage	BEAUTY		P 36	
Abdominal area drainage	BEAUTY		P 37	
Upper limb drainage	BEAUTY		P 38	
Leg firming	BEAUTY		P 39	P15
Hip-gluteus firming	BEAUTY		P 40	
Abdominal area firming	BEAUTY		P 41	

Upper limb firming	BEAUTY		P 42	
Leg toning	BEAUTY		P 43	P 16
Dorsal toning	BEAUTY		P 44	
Hip-gluteus toning	BEAUTY		P 45	
Abdominal area toning	BEAUTY		P 46	
Upper limb toning	BEAUTY		P 47	
Leg sculpting	BEAUTY		P 48	P 17
Hip-gluteus sculpting	BEAUTY		P 49	
Abdominal area sculpting	BEAUTY		P 50	
Pectoral sculpting	BEAUTY		P 51	
Sculpting	BEAUTY	P 4		
Leg vascularization	BEAUTY		P 52	
Upper limb vascularization	BEAUTY		P 53	
Leg lipostress	BEAUTY		P 54	
Gluteus lipostress	BEAUTY		P 55	
Hip lipostress	BEAUTY		P 56	
Abdominal area lipostress	BEAUTY		P 57	
Upper limb lipostress	BEAUTY		P 58	
Leg mass building	BEAUTY		P 59	
Upper limb mass building	BEAUTY		P 60	
Abdom. area post-pregnancy drainage	BEAUTY		P 61	P 18
Abdom. area post-pregnancy lipolysis	BEAUTY		P 62	P 19
Abdominal area post-pregnancy toning	BEAUTY		P 63	P 20
Man abdominal area definition	BEAUTY		P 64	
Man pectoral definition	BEAUTY		P 65	
Woman leg mass building	BEAUTY		P 66	
Woman upper limb mass building	BEAUTY		P 67	
Swollen upper limbs	BEAUTY		P 68	
Swollen legs	BEAUTY		P 69	

Leg toning massage	BEAUTY		P 70	
Face beauty 1	BEAUTY		P 71	
Face beauty 2	BEAUTY		P 72	
Breast toning 1	BEAUTY		P 73	
Breast toning 2	BEAUTY		P 74	
Conventional antalgic Tens	TENS	P 5	P 75	P 21
Endorphinic Tens	TENS	P 6	P 76	P 22
Rotator cuff tendinitis	TENS	P 7	P 77	P 23
Muscle pain	TENS	P 8	P 78	P 24
Knee osteoarthritis	TENS	P 9	P 79	P 25
Post-surgery pain	TENS	P 10	P 80	P 26
Muscle injury	TENS	P 11	P 81	P 27
Cervical pain	TENS	P 12	P 82	P 28
Trapezius pain	TENS	P 13	P 83	P 29
Sciatica	TENS	P 14	P 84	P 30
Lumbago	TENS	P 15	P 85	P 31
Epicondylitis	TENS	P 16	P 86	P 32
Scapulo-humeral periarthritis (shoulder)	TENS	P 17	P 87	P 33
Carpal tunnel	TENS	P 18	P 88	P 34
Osteoarthritis	TENS	P 19	P 89	P 35
Menstrual pain	TENS	P 20	P 90	P 36
Ankle prevention	PREVENTION		P 91	P 37
Knee prevention	PREVENTION		P 92	P 38
Quadriceps muscle prevention	PREVENTION		P 93	P 39
Lower limb reathletization	PREVENTION			P 40
Shoulder reathletization	PREVENTION		P 94	P 41
Shoulder prevention	PREVENTION		P 95	P 42
Elbow prevention	PREVENTION		P 96	P 43

Upper limb reathletization	PREVENTION			P 44
Ankle muscle tone recovery	PREVENTION			P 45
Vastus medialis reinforcement	REHAB			P 46
Swollen ankles	REHAB			P 47
Hemiplegia upper limbs	REHAB			P 48
Hemiplegia lower limbs	REHAB			P 49
Recovery after ACL surgery	REHAB			P 50
Functional recovery lower limbs	REHAB			P 51
Quadriceps atrophy with prosthesis	REHAB		P 97	P 52
Shoulder subluxation prevention	REHAB		P 98	P 53
Muscle spasms	REHAB			P 54
Multiple sclerosis - Leg flexor extensor mm.	REHAB			P 55
Multiple sclerosis - Ankle flexor extensor mm.	REHAB		P 99	P 56
Mixed incontinence	INCONTINENCE		P 100	P 57
Stress incontinence	INCONTINENCE			P 58
Urge incontinence	INCONTINENCE			P 59
Ionophoresis	IONO			P 60
	TOTAL	20	100	60
	SPORT	3	30	12
	BEAUTY	1	44	8
	TENS	16	16	16
	PREVENTION	0	6	9
	REHAB	0	3	11
	INCONTINENCE	0	1	3
	IONO	0	0	1
	Total	20	100	60

**Because of the presence of clinical programs, this product is a medical device. Therefore*

it is certified by the KIWA CERMET ITALIA n°. 0476 according to 93/42/CEE directive for medical devices. The certification covers clinical applications.

GENERAL NOTES ON ELECTRODE POSITIONING

Correct electrode positioning and size choice are fundamental to assure the effectiveness of electrostimulation.

To choose the size of the electrodes and their positioning it is necessary to refer to the images at the end of this manual. Information is also available on our website www.globuscorporation.com.

N.B. For all the programs that cause significant muscle contraction (e.g. strength, hypertrophy, toning and firming programs) it is fundamental to place the electrode on the muscle **motor point**, which is the most sensitive to stimulation.

If the electrode is not positioned exactly on the motor point, the contraction could be small and/or annoying. In this case it is necessary to shift the positive electrode of a few millimeters to feel an effective and comfortable muscle contraction.

Body position during stimulation


Body position during electrostimulation depends on the body part involved and on the type of program that is being carried out.

During high-intensity treatments, we suggest blocking the limbs in order to work isometrically. For instance, if you want to treat quadriceps with a strength program, we suggest to carry out the treatment in a sit position and block the feet, in order to avoid involuntary leg extension during the contraction.

For all the programs that do not imply high intensity (massages, decontracting, drainage...) the body position is not important, as long as it is comfortable.

Electrode positioning for Tens programs

In the following pages of this manual you can find the images with the correct electrode positioning for tens treatments. If the localization of your pain type is not included among the images, you can place the electrodes by forming a "square" on the aching area. Here you have an example.

TENS

(use colored cables)

Indications for the use of ionophoresis

This program uses a low-intensity and continuous current (between 5 and 10 mA) that helps the absorption of the active ingredients of a drug in a specific area of the body. It is often required that a drug is administered directly on the area that needs to be healed, avoiding an oral or intravenous administration. With ionophoresis it is possible to administrate a drug through transcutaneous penetration directly on the body area that needs to be treated.

The drug must be prepared in ionic form and must have an electric charge. The principle on which ionophoresis is that the continuous current flows through skin from one electrode to the other (from anode to cathode) consequentially transporting the electrically charged ions that are inside the treated body area.

In practice, sponges soaked in a solution containing the active ingredient are used, bearing in mind that if it has a negative charge, it must be applied to the negative electrode and vice versa.

Please note that, in the cables supplied, the white output corresponds to the positive pole and the red to the negative pole and to use only channel 1 of the device. The current will carry the active principle inside the tissues, because the ions of the active principle itself will migrate towards the opposite pole until the product is completely absorbed. The kit with sponges and electrodes for ionophoresis is an extra accessory that can be purchased separately.

The intensity value is adjustable from a minimum of 0 mA to a maximum of 10 mA. To ensure safety, the maximum intensity, calculated with the smallest ionophoresis electrodes (50 x 50 mm), is 0.40 mA per cm² of electrode.

DENSITY = Current displayed in mA / Electrode area in cm²

Indications for the application:

- The manufacturer declines the responsibility of the prescribing physician for the choice of the drug to be used with ionophoresis. For this purpose, in this user manual the manufacturer provides the data for a correct definition of the protocol and the drug that has to be used.
- Before using ionophoresis, therefore, it is mandatory to contact the prescribing doctor.
- Both sponges in the ionophoresis kit must be wetted with distilled water or physiological solution, then insert the silicone electrodes inside the sponges and connect the cables.
- The drug used for the therapy should never be applied directly on the skin but to the absorbent surface of the electrode corresponding to the polarity of the drug itself. In fact, it is very important to consider the positive or negative charge sign of the active ingredient of the drug, which must be applied correctly: the drug must be placed on the electrode with the same charge sign as the active ingredient.

- Place the two electrodes on the area that must be treated, at a distance of about 10-20 cm and fix them with the elastic band.
- Gradually increase the intensity until the patient feels a slight tingling.
- Once the program is finished, turn the device off, and disconnect the electrodes. Clean the sponges and bands following the instructions on the drug leaflet. If not indicated, it is however recommended to wash the sponges and elastic bands very well, with plenty of hot water and soap, so that the next time there are no traces of the drug.

Indications for the treatment with incontinence programs

The use of muscle stimulation programs referred to the treatment of urinary incontinence can be used exclusively for the treatment of urge, stress and mixed incontinence.

These three programs are indicated for the reinforcement of the pelvic floor, that often cause incontinence. The program must be chosen after a diagnostic of the type of incontinence (urge, stress and mixed). Urological applications entail the use of vaginal or anal probes for the specific intended use, which must be covered by the CE MDD certification, in compliance with the directive 93/42/CEE. Such probes must be bipolar and equipped with a 2 mm female connector for 2 mm male cables. Such probes are available in stores with the following reference codes:

REF G1188 Single-patient vaginal probe


REF G0757 Single-patient anal probe.

Warnings

Since these applications have a medical intended use, they must be carried out upon previous medical authorization.

Use

- To use the probes correctly, follow the instructions provided by the manufacturer or the physician.
- Then connect the stimulation cable to the probe: the connector of the stimulation cable with the + symbol must be connected to the red connector of the probe. The - symbol must be connected to the black connector of the probe.

- Connect the cable to the device.
 - For the vaginal stimulation, lie flat supine with bended knee and opened legs. To use the anal probe, lie on one side in lateral decubitus position
 - Insert the probe in the vagina or anus using, if necessary, some water-based gel
 - Choose and start the incontinence program
 - Gradually increase the intensity until perceiving a muscular contraction.
- Once the program is finished, turn the device off, extract the probe and disconnect the connection cables. Clean the accessory and put it in a clean bag.

TREATMENT DESCRIPTIONS

Following is a short description of the programs provided by the electrostimulator. For further information about the programs you can download from our website a complete guide, where you can find all the indications and suggestions to carry out the treatment in the correct way.

SPORT

WARM-UP: Increases muscular temperature for better training preparation. Increases the blood flow to and metabolism of the muscle. Recommended before training sessions or competitions.

Intensity: low or intermediate.

MAXIMUM STRENGTH: Improves maximal force capability and increases muscular mass. Allows training with maximum load and reduced trauma risk as compared to classical training. Recommended for activities requiring great physical strength. Also useful during the early preparation period as it trains the muscles to receive significant stress.

Intensity: maximum, at pain threshold limit. During active Rest, adjust the intensity to obtain comfortable contractions.

ENDURANCE STRENGTH: Improves the ability to display a high level of force for an extended period of time and to resist toxin accumulation, thus delaying muscular fatigue. The program involves a long series of lengthy contractions, alternated with brief active rest. Recommended for sports requiring intense muscular work for an extended period of time.

Intensity: intermediate or maximum. Active Rest should induce muscular vibrations.

EXPLOSIVE STRENGTH: The Explosive Force program increases the ability to rapidly achieve maximum levels of force. The objective is to recruit the maximum number of muscle fibers in the shortest time possible. It is recommended to begin with a Maximal Force program to prepare the muscle for the Explosive Force program. Recommended in all sports requiring maximal force under natural load conditions,

such as sprints, ball sports, and alpine skiing. Intensity: maximum muscular contraction at pain threshold limit. Active Rest should induce muscular vibrations.

DECONTRACTING: Enables effective muscle relaxation thanks to the pleasant contractions that increase blood circulation. These programs can be used whenever muscular relaxation is needed. Intensity: low to intermediate.

CAPILLARIZATION: Increases the blood flow to the muscle and improves muscle resistance and recovery. An increased arterial capacity keeps the primary and secondary capillary network active, thus improving the tissue oxygenation system. This delays fatigue during very intense work. Capillarization programs are particularly recommended for sports requiring resistance force and endurance.

Intensity: intermediate.

POST-COMPETITION/TRAINING ACTIVE RECOVERY: Improves recovery from muscular fatigue after training or competitions and enables easier training resumption. The stimulation varies from muscular vibrations to light contractions. This allows greater circulation and reduced muscular hardening. It also alleviates pain thanks to the release of endorphins. This program is valid for all sports as it improves and speeds up the recovery of maximum muscular capacity. This program is recommended in the hours immediately after training or competitions.

Intensity: low to intermediate.

BEAUTY

DRAINAGE: The coordinated stimulation obtained by using all channels of the electrostimulator produces an in-depth muscular massage effect. The massage activates venous and lymphatic systems by means of a mechanical effect. The purpose is to improve blood circulation and enable the toxin elimination. Intensity: low to intermediate.

LIPOLYSIS: The low frequency of this program positively affects the metabolism and increases the elimination of toxins and excess fluid. This is a method for reshaping critical areas, modifying the metabolism of lipocytes. Improves blood circulation. Helps to eliminate the orange peel skin. Intensity: low to intermediate.

VASCULARIZATION: Enables effective muscle relaxation thanks to the pleasant contractions that increase blood circulation. These programs can be used whenever muscular relaxation is needed. Intensity: intermediate.

FIRMING: This program is suitable for sedentary individuals with flaccid muscles that require reactivation of muscular tone. It prepares the muscles for more intense stimulations. It can be performed every day on large muscle groups such as thighs, abdominal and dorsal muscles. Intensity: intermediate to high.

TONING: The toning program allows to develop and maintain muscular tone. This type of treatment may be combined with aerobic activity and is suitable for users of any level. The program is especially directed at all those who, mindful of their

own physical condition, wish to integrate their own regular training. The program can be used up to three times a week. Intensity: intermediate to high.

SCULPTING: The sculpting program is a merely aesthetic treatment and is ideal for those who wish to better define their already toned muscle system. The work performed is very intense, so slight muscle soreness may be experienced the day after the treatment. However, this effect may be avoided by performing stretching exercises before and after the session or a winding- down program. The program is suitable for those who wish to achieve muscular definition in those areas where traditional muscular development exercises are insufficient. Intensity: intermediate to high.

MASS BUILDING: A specific program for those who wish to increase muscle mass and volume. Intensity: maximum muscular contraction at pain threshold limit.

A muscle warm-up program should be carried out before performing any of the following programs: Firming, Toning, Sculpting and Mass Development.

TENS

Uses mechanisms of the nervous system to alleviate both acute and chronic pain. Impulses travel along the nerves to block the pain propagation paths. This analgesic effect occurs during program execution and can last for a few hours. However, its duration may vary. The TENS programs also improves circulation. Place the electrodes on the painful area. These programs can be used in any situation requiring pain relief. Caution! TENS programs alleviate the pain but do not eliminate the cause. If the pain persists, seek medical advice.

TENS (Antalgic-Endorphinic): Allows pain relief and effective muscle relaxation. TENS programs use mechanisms of the nervous system to alleviate both acute and chronic pain. Impulses travel along the nerves to block the pain propagation paths. Muscular relaxation improves blood circulation. It can be used in any situation requiring both pain relief and muscle relaxation. Caution! TENS programs alleviate the pain but do not eliminate the cause. If the pain persists, seek medical advice. Intensity: ANALGIC TENS: tingling but not painful sensation. Should not induce muscular contractions. ENDORPHINIC TENS: visible muscular vibrations.

TENS (Cervical pains, Lumbago, Periarthritis, Epicondylitis, Sciatica, ...): Uses mechanisms of the nervous system to alleviate both acute and chronic pain. Impulses travel along the nerves to block the pain propagation paths. This program combines stimulation with light contractions that stimulate the release of endorphins. Caution! TENS programs alleviate the pain but do not eliminate the cause. If the pain persists, seek medical advice. Intensity: tingling but not painful sensation.

REHAB

VASTUS MEDIALIS REINFORCEMENT AND RECOVERY AFTER ACL SURGERY: Resumption of physical activity after a period of interruption. After a trauma or other interruptions of physical activity, muscles may lose their tone and mass. In these cases, normal muscular capacity must be regained as soon as possible. The program enables training of specific muscles concerned.

It can be used in any situation requiring activity resumption. Intensity: low to intermediate.

PREVENTION PROGRAMS: The stimulation is aimed at neuromuscular facilitation, which especially improves the mobilization of the injured articulation. Prepares the limb to resume its normal functionality. It is recommended to use this program during the first phases of the recovery of muscular trophism. Intensity: intermediate.

RE-ATHLETIZATION: Increases the diameter and capacity of muscle fibers.

This program is recommended after using the Atrophy program for a period, as soon as the muscle shows a slight recovery of volume and tone. Intensity: Intermediate.

HEMIPLEGIA: If the patient suffers from a perception problem with decreased attention (semi-inattention), integrating the treatment with training of this ability may be extremely useful. The treatment must not be merely passive, but the patient can take on an active role in the training. The program can be used for the facilitation and rehabilitation of motor faculties as well as a method to reduce spasticity.

Intensity: depending on the muscle district to be stimulated.

INCONTINENCE: Programs recommended for the stimulation of perianal muscles.

These programs are recommended in case of urological pain or incontinence.

Programs to be used with probes in compliance with the law.

ACTION PRINCIPLES

Muscular electrostimulation

Electrostimulation is a technique that, by means of electric pulses that act on the muscle motor points (motoneurons), causes a muscular contraction similar to voluntary contraction.

The majority of human muscles are striated or voluntary, with approximately 200 muscles on each side of the body (about 400 on the whole).

Electrostimulation advantages

The electro stimulation does not want to replace physical activity but it must be considered as an integrating treatment.

According to the aims it will be useful:

- To the sportsmen/sportswomen, in order to integrate the normal training and to increase performances
- To who wants to intervene on some imperfections
- To who suffers from some pathologies
- To who has suffered from a trauma or in the rehabilitation

The electro stimulation is a technique that provokes a muscular contraction completely similar to the volunteer one, using electrical TENS impulse. There are two different ways of use:

- Muscular stimulation (ideal for the development of the force and for aesthetic treatments)
- The stimulation to nervous ends (ideal for the treatments against the pain).

Types of muscle

The muscle can be subdivided in three different types: streaked or voluntary muscle; cardiac muscle and smooth or involuntary muscle.

The voluntary muscle includes the skeletal muscle that gives movement to the skeleton and the cutaneous muscle that acts on facial expressions.

The cardiac muscle and the smooth muscle are not voluntarily controlled. The most part of the muscles of the human body belongs to the category of streaked or voluntary muscles, with approximately 200 muscles for every side of the body (approximately 400 in total).

Skeletal muscles are the target of the EMS.

Mechanism of the muscular contraction

The skeletal muscle exercises its functions through the mechanism of the contraction. When the muscular contraction happens, the movement of the articulations is produced and consequently the movement of the skeleton, too.

The muscle contracted in the following way: when a person decides to make a movement, an alteration of the electric potential in the motor center of the brain and an electric impulse is transmitted to the muscle that has to move.

The electrical impulse is transmitted as a change of potential (voltage); an electrical impulse crosses the motor nerve to the muscle that must be stimulated. After the reception of the impulse, Ca ions are freed by the reservoirs causing the approach of actin and myosin. The result is that the distance between the bands Z diminishes causing therefore the muscular contraction.

The energy demanded for the contraction is supplied by the supply of sugars and fats stored in the human body. In other words, the electrical stimulation is not a directed energy resource but it works as an instrument that triggers the muscular contraction. The same type of mechanism is activated when the muscular contraction is produced by the EMS.

They assume in other words the same role of a natural impulse transmitted by the motor nervous system. Normally the muscle relaxes and returns to its original state them at the end of the contraction.

Isotonic and isometric contractions

The isotonic contraction reveals when, in a movement of our body, the interested muscles produce a state of constant tension. The isotonic contraction reveals when, in a movement of our body, the interested muscles produce a state of constant tension. Instead, when the joint heads of a physical segment are blocked and the musculature produces a tension, we talk about an isometric contraction. In the event of the electro stimulation an isometric contraction is preferred because it fosters to obtain a more effective control.

The distribution of the different types of fibers in the muscle

The relationship between the two main categories (type I and the type II) can vary in sensitive way. There are muscular groups that are typically constituted by fibers of type I, like the soleus, and there are muscles that have only fibers of type II like the orbicular muscle. But in the most part of the cases we have a simultaneous presence of various types of fibers.

The studies lead on the distribution of fibers in the muscle have put in evidence the strait relationship that exists between the motoneuron (tonic or phasic) and the functional features of fibers from it innervate and have demonstrated as to a specific motor activity (and sports in particular) can determine a functional adaptation of the fibers and a modification of the metabolic features of the same ones.

Motor unit type	Contraction type	Contraction frequency
Tonic ST	slow contraction I	0 - 50 Hz
Phasic FT	fast contraction II	50 - 70 Hz
Phasic FTb	fast contraction II b	80 - 120 Hz

In order to make a tissue pass from the phase of rest to the one of excitation, by means of an induced electrical stimulus (impulse of the electrostimulation), some conditions are necessary:

- the current amount
- the duration of the stimulus that must be adapted to the corporeal district that is wanted to be stimulated.

From this consideration, for the excitation of a tissue, the relationship between the duration and the intensity of the stimulus is particularly important and not only the value of the peak of the intensity. This relationship varies according to the muscular district.

Stimulation intensity

The current intensity necessary to obtain muscular contraction is personal and depends on the position of the electrodes, the underlying adipose tissue layer, sweating, the presence of hair on the area to be treated, etc.. For these reasons, the same current intensity can generate different feelings from person to person, from day to day, and from the right side to the left side of the body. During different sessions, it will be necessary to regulate the intensity in order to obtain the same level of contraction because of the accommodation phenomena. The current intensities recommended in the different phases are proposed as indicative values, and each person should modify these levels according to his/her personal needs.

- Moderate: the muscle does not tire, even during long treatment. The contraction is agreeable and tolerable. First level of the intensity graph.
- Intermediate: the muscle is visibly contracted but the stimulation does not trigger the joint movement. Second level of the intensity graph.
- High: the muscle is contracted noticeably. The muscle contraction would extend or bend the limb if not blocked. Third level in the intensity graph.
- Maximal: the muscle is contracted maximally. This is an intense treatment that should be performed only after many applications.

Moderate	From 10 mA to 20 mA
Intermediate	From 20 mA to 30 mA
High	Above 30 mA
Maximal	On the verge of the tolerance limit, always under the threshold of pain

In the descriptions of the treatments, the best levels of intensities are recommended. NOTE: The recommended levels of current are only indicative.

Tens

Transcutaneous Electrical Nerve Stimulation (TENS) is a selective stimulation of the large fibers of the peripheral nerves favoring the closing of the gate entrance for the pain pulses and increasing the release of endorphin substances, reducing in this way the pain intensity. Therefore with TENS we want to treat the severe and chronic pain due to the main musculoskeletal pains.

The pain decrease following to the TENS currents application is due to these factors:

- a. Gate control theory
- b. Endorphin secretion
- c. Different sedative effects in relation with the frequency

Gate theory

If the electric signals that lead to the brain information about pain are stopped, also the pain perception is eliminated. If, for instance, we hit our head into an object the first thing we do is massaging the area affected by the trauma. In this way we stimulate the receptors related to touch and pressure. TENS in continuous mode and in frequency modulation can be used to generate signals similar to the ones of touch and pressure. If their intensity is enough, their priority is so high that it prevails on the pain signals. Once the priority is gained, the gate related to the sensory signals is opened and the pain one is closed, impeding in this way the passage of these signals to the brain.

Endorphin secretion

When a nervous signal proceeds from the pain area to the brain, it spreads through a chain of connections joined together called synapse. The synapse can be seen as the space between the end of a nerve and the start of another. When an electric signal reaches the end of a nerve, it produces some substances called neurotransmitters that pass through the synapse and activate the start of the next nerve. This process repeats for all the length needed to the signal to reach the brain. The opioids involved in the pain reduction have the task to insinuate in the synapse space and impede the neurotransmitter propagation. In this way a chemical block of the pain signals occurs. The endorphins are opioids naturally produced by the body to tackle the pain and they can act both on the marrow and on the brain, in this way they are effective analgesics. Tens can increase the natural production of endorphins and, thereby, they act decreasing the pain perception.

Incontinence programs

The urology programs require the use of specific endovaginal and endorectal electrode probes, certified according to the Directive for Medical Devices 93/42/EEC. These are bipolar probes with a 2-mm female adapter which attaches to 2-mm male cables.

Warnings

Urological electrostimulation is a medical application, which must be carried out under medical supervision.

Use

To correctly use the probe electrode, follow the instructions provided by the manufacturer and given by the physician.

Maintenance

For cleaning, sterilization and disinfection, refer to the manufacturer's instructions.

Suggestions

In case of deterioration of the probe electrode, replace it immediately.

Ionophoresis

Ionophoresis is an electrotherapy form that acts transmitting pharmacological substances inside the tissue thanks to a continuous electric current of unidirectional type.

Ionophoresis bases on the capacity of ionic dissociation of some medicated substances, of very low molecular weight, after being dissolved in water.

It is really important to know that the active part of the medicine has, after being dissociated under ionic form, positive or negative charge, with the aim to place it correctly according to the direction of the electric flux.

The ions of the medicated substance are transmitted inside the organism through cutaneous areas that oppose a low resistance to the current reaching in this way the cellular membranes that are thereby electrically modified.

The home user can use the ionophoresis treatments only after consulting the specialist who will prescribe the medications to use and give the indications for the type of currents to use.

The current intensity should be regulated as to be barely perceptible. We recommend starting the treatment with the lowest intensity and gradually increasing.

DO NOT APPLY THE MEDICATION DIRECTLY TO THE SKIN. Apply the medication to the absorbent surface of the electrode corresponding to the medication's polarity; the absorbent surface of the other electrode should be dampened with slightly salted water, to promote conductivity.

MAINTENANCE AND CLEANING

Maintenance and cleaning of the device

- In case of real or alleged malfunction, do not tamper with the device or try to repair it yourself.

Do not intervene on the device and do not open it. Only specialized and authorized centers can repair it.

- Avoid violent impacts that may cause damage and malfunctions to the device even if undetectable immediately.

- Use this device in a dry environment and in an open space (not wrapped in any materials).

- Clean the device and accessories only with disinfectant with sodium hypochlorite or quaternary ammonium salt diluted with distilled water equal to 0,2-0,3%. After cleaning/disinfecting it, dry the device and its accessories with a clean cloth.

- Always use the device and its accessories with clean hands.

- It is recommended to use the device in a clean room, to avoid the contamination of the device with dust and dirt.
- It is recommended to use the device in a ventilated space, with regular air change.
- It is recommended to clean/disinfect the parts after every use.

Accessories

After using the multi-purpose single patient and/or single-use electrodes, they must be stored using their plastic film and placed in the plastic bag.

Avoid that the electrodes touch each other or that they overlies one over the other. Once the package has been opened, the electrodes can be used for 25-30 applications. The electrode must be always used with clean hands and they must be replaced if they are not perfectly in contact with the skin.

If using non self-adhesive electrodes it is suggested to clean the surface with proper cleansers that respect the requirements described in the manual.

The electrodes must be stored in their bag and in an environment that respects the requirements described in the manual.

After the end of a treatment, unplug the cables from the connectors and clean them carefully with proper cleansers that respect the requirements described in the manual. After cleaning and drying them, they must be folded up and placed in the plastic bags supplied with the cables.

Battery: how to charge the batteries

The device is supplied with a set of rechargeable nickel-metal hydrate batteries (4,8V 800mAh).

When in the display the battery to be charged icon () appears, we suggest recharging the device, as after the appearance of this symbol the electrostimulator autonomy is limited and could not allow the execution of a complete program.

To charge the batteries, turn off the electrostimulator and disconnect the electrodes, then connect the electrostimulator to the charger provided by inserting the plug in the appropriate inlet.

We suggest charging the device for 8-12 hours to obtain the greater duration of the battery pack.

The device will be operating for 15-20 hours after the charge, according to the treatments carried out. It is better to recharge the device during this time, which corresponds to 70 treatments approximately.

A possible increase in the temperature in the battery area during the charge is normal. Do not use a different charger from the one provided with the device. To replace the electrostimulator batteries, contact an authorized service center.

WARRANTY

The device is guaranteed to the first user for a period of twenty-four (24) months from the date of purchase against defects in materials or of the manufacturing (if the user uses the device for professional purposes, is twelve, 12, months), provided that it is used properly and maintained under normal operating conditions.

Warranty coverage is limited in the following cases:

- Six (6) months for supplied accessories such as batteries, chargers, power supplies, cables, etc.;
- No warranty for extendable accessories and materials such as, for example, electrodes, etc..

The warranty is valid and enforceable in the country where the product was purchased. In the event that the product is purchased in any country of the European Community, the warranty is still valid in all its countries.

In order to take advantage of the warranty service, the user must comply with the following warranty clauses:

1. The products, and all accessories, to be repaired must be sent by and at the expenses of the customer in their original packages.
2. The product's warranty is subject to the production of a fiscal document (fiscal receipt, receipted bill or sale invoice), attesting the product's purchase date.
3. The repair work shall have no effect on the original warranty expiry date and shall neither renew nor extend it.
4. If, during the repair work, no defects are found, the costs related to inspection times shall in any case be charged.
5. The warranty becomes void where the fault has been caused by: impacts, falls, erroneous or improper use of the product, use of non-original power supply unit or charger, accidental events, alteration, replacement/detachment of the warranty seals and/or tampering with the product. The warranty does not cover damages caused during transportation when unsuitable packages are used.
6. The warranty does not cover the inability to use the product, other incidental or consequent costs or other expenses incurred by the purchaser.

The features and dimensions reported in this manual are not binding.

The manufacturer reserves the right to make changes without previous notice.

FREQUENTLY ASKED QUESTIONS

Why, when I increase the current intensity above 1 mA, it returns immediately to 0 mA?

Before starting any program, check that the cables are connected and that the electrodes are placed on the area to be treated. In addition, be sure that the electrodes are not worn, because this could inhibit their ability to conduct current.

What kind of electrodes should be used for electrostimulation?

Use self-adhesive electrodes, which are practical and improve the quality of stimulation. If used with care, they will last for 25-30 applications. The electrodes should be replaced when they have no longer a good connection with the skin.

Where should the electrodes be placed?

In the back of this manual, there is a comprehensive electrode placement guide (it is not necessary to respect the polarities indicated). You can follow this instructions. To verify the correct placement of the electrodes, use the special Find Motor Point Pen program or follow this empirical method: place the electrodes as indicated in the pictures in the back of this manual; start the stimulation; with your hand, move the electrode by sliding it along the muscle without removing it from the skin. You will notice a change in contraction according to the different positions. Once you locate the point where the stimulation is greatest, decrease the channel intensity to zero (0,0 mA), replace the electrode and increase gradually the intensity.

Do the splitter cables and the split cables allow you to use more electrodes with the same channel?

This permits, for instance, to stimulate vastus medialis and vastus lateralis of the quadriceps with one single channel. They are not recommended for medical applications.

Does the power decrease using Y cables?

The power intensity for each channel does not vary. However, when Y cables are used to split one single channel in two, the current is distributed on a wider muscle area, therefore contraction will be less pronounced. Increase the intensity to obtain the same contraction level.

Can electrostimulation hurt me?

It is very unlikely that electrostimulation damages muscles. In this case it is important to increase intensity gradually while observing the muscle behavior and avoiding to keep the limb completely outstretched. In case of doubt, please contact a specialist.

Is it possible to use the electrostimulator during menstruation cycle?

Electrostimulation may interfere in some way with menstruation, causing anticipation, delay, accentuation or reduction of the cycle; however, these effects

are subjective and highly variable. It is recommended to avoid treatments in the abdominal zone during menstruation cycle and immediately before or after it.

Is it possible to use the electrostimulator during lactation?

Until now, no collateral effects regarding lactation have been observed. Nonetheless, during lactation, it is recommended not to stimulate the thoracic region.

Are dermatological diseases (e.g. psoriasis, urticaria) contraindications for electrostimulation?

Yes, do not treat areas affected by dermatological diseases.

When are the first results visible?

Aesthetic results of electrostimulation are always subjective. For Toning program, with a regular program of 3-4 sessions per week, a noticeable result may be observed after 15 days. For Lipolysis and Drainage programs, 40 days of treatment are necessary. Results are obtained more quickly if treatments are combined with good physical activity and a correct life style.

How many sessions can I weekly perform?

For physical training, consult the program of weekly training in the Globus Personal Trainer. For fitness and aesthetics programs, the number of sessions depends on the type of treatment: 3-4 sessions per week on alternate days are suggested for toning; daily treatments are permitted for Lipolysis and Drainage programs.

EMC accompanying documents

Essential performance

PERFORMANCE	CONDITION	RISK	ACCEPTED EVENT
Electrostimulation.	External disturbance (Burst).	Display information no longer readable.	The machine must stop the stimulation.
			The machine must maintain the stimulation and accept the commands.
	Lack of internal power supply.	Interruption of the treatment.	The machine must signal the battery exhaustion and the interruption of the treatment.
	Lack of external power supply.	Interruption of the treatment.	The device, if equipped with a battery, must continue the treatment signaling that operation is carry out in battery mode.
	Detachment of an electrode.	Unpleasant stimulation or painful electric shock in case of	The device must constantly monitor the current on each active channel set over 9 mA. In case the detected

		reconnection of the electrode.	current is below a certain threshold, the machine must rest the current of the channel.
	Setting of a current that is too high	Dangerous stimulation.	The device must derate the voltage boost stage to prevent emitting a current beyond the maximum value. The machine that does not detect the electrode contact must reset the channel current and must indicate this reset in the indication on the display.
Loading of the programs from the memory.	Error in the data from the memory.	Execution of an incorrect program.	The machine must check the correctness of the data of the programs. In case an error is detected, the device must restart.
Battery charge.	Battery overheating.	Damaging of the device, display information no longer readable, explosion, fire.	The device must monitor the temperature of the battery, if a certain threshold is exceeded, the battery charging must be interrupted.

In compliance with:

EN 60601-1: 2006 + A1: 2013

EN 60601-1-2: 2015

EN 60601-2-10: 2015

EN 60601-1-11: 2015

Warning: radiofrequency communication devices (including accessories like antennas or antenna cables) must be used at least 3 meters away from every part (including cables and accessories) of the device. Otherwise performance can be affected.

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