

LIPOZERO G150

AESTHETIC CAVITATION | 42 PROGRAMS

ULTRASONIC CAVITATION IS THE MOST EFFICIENT METHODOLOGY TO FIGHT CELLULITE AND TO ELIMINATE EXCESS FAT. THANKS TO LIPOZERO G150 YOU CAN BENEFIT OF THE PROFESSIONAL TREATMENTS OF AN ESTHETIC CENTER COMFORTABLE AT HOME.

10 MF-PLUS | 10 CAVIMAX 150kHz | 8 BEAUTY | 8 SKIN | 5 ADIPOSE MOBILITATION NorAd-526 | 2 MAN ADIPOSITIVITY

TECHNICAL FEATURES

Display

Visible area size 2,6"

Output

1

Head

1 head Ø mm 42 - IP68

Power

3W/cm² ± 20%

Emission Frequency

3MHz - 1MHz - 526KHz - 150KHz

Emission

Continuous and pulsed

Power supply

- Mains
- Rechargeable batteries (not included)

Size


mm 170x220x60

Weight

gr 870

EQUIPMENT

- 1 Carrying case
- 1 Ultrasound LIPOZERO G150
- 1 Multi-frequency head
- 1 Gel flask
- 1 Power supply unit
- 1 Operating manual



CE



MAINS SUPPLY



RECHARGEABLE BATTERIES

GLOBUS
ITALIAN EXCELLENCE

LIPOZERO G150

AESTHETIC CAVITATION | 42 PROGRAMS

PROGRAMS LIST

CAVIMAX 150KHZ

Lipolysis

Abdominal lipolysis

Cellulite

Oedematous cellulite

Fibrous cellulite

Culotte de cheval

Cellulite oedematous legs

Cellulite fibroses legs

Cellulite with orange skin effects

MF-PLUS

Lipolysis

Cellulite

Cellulite oedematous

Cellulite fibroses

Culotte de cheval

Cellulite with orange skin effects

Man adiposity abdomen

Man adiposity hip

Man adiposity lower limbs

Man adiposity upper limbs

ADIPOSE MOBILISATION NorAd-526

Lipolysis - NorAd-526

Cellulite - NorAd-526

Adiposity - NorAd-526

NorAd-526 1

NorAd-526 2

MAN ADIPOSITY

Man adiposity abdomen

Man adiposity hip

BEAUTY

Massage

Modelling

Modelling legs

Drainage

Drainage legs

Drainage arms

Stretch marks

Active principal application

SKIN

Acne

Face lifting

Wrinkle treatment

Scar

Couperose

Body spots

Face spots

Restructuring skin