

# LIPOZERO G150

## AESTHETIC ULTRASOUND | 42 PROGRAMS

ULTRASONIC CAVITATION IS THE MOST EFFICIENT METHODOLOGY TO FIGHT CELLULITE AND TO ELIMINATE EXCESS FAT. THANKS TO LIPOZERO G150 YOU CAN BENEFIT OF THE PROFESSIONAL TREATMENTS OF AN ESTHETIC CENTER COMFORTABLE AT HOME.

10 MF-PLUS | 9 CAVIMAX 150 kHz | 8 BEAUTY | 8 SKIN | 5 ADIPOSE MOBILITATION NorAd-526 | 2 MAN ADIPOSITIVITY

### TECHNICAL FEATURES

#### Display

Visible area size 2,6"

#### Output

1

#### Head

1 head Ø mm 42 - IP68

#### Power

3 W/cm<sup>2</sup> ± 20%

#### Emission Frequency

3 MHz - 1 MHz - 526 kHz - 150 kHz

#### Emission

Continuous and pulsed

#### Power supply

- Mains
- Rechargeable batteries (not included)

#### Size


170x220x60 mm

#### Weight

870 gr

### EQUIPMENT

- 1 Carrying case
- 1 Ultrasound LIPOZERO G150
- 1 Multi-frequency head
- 1 Gel bottle
- 1 Power supply unit
- 1 Operating manual



CE



MAINS SUPPLY



RECHARGEABLE BATTERIES

**GLOBUS**  
ITALIAN EXCELLENCE

# LIPOZERO G150

## AESTHETIC ULTRASOUND | 42 PROGRAMS

### PROGRAMS LIST

#### CAVIMAX 150KHZ

Lipolysis

Abdominal lipolysis

Cellulite

Oedematous cellulite

Fibrous cellulite

Culotte de cheval

Cellulite oedematous legs

Cellulite fibroses legs

Cellulite with orange skin effects

#### MF-PLUS

Lipolysis

Cellulite

Cellulite oedematous

Cellulite fibroses

Culotte de cheval

Cellulite with orange skin effects

Man adiposity abdomen

Man adiposity hip

Man adiposity lower limbs

Man adiposity upper limbs

#### ADIPOSE MOBILISATION NorAd-526

Lipolysis - NorAd-526

Cellulite - NorAd-526

Adiposity - NorAd-526

NorAd-526 1

NorAd-526 2

#### MAN ADIPOSITY

Man adiposity abdomen

Man adiposity hip

#### BEAUTY

Massage

Modelling

Modelling legs

Drainage

Drainage legs

Drainage arms

Stretch marks

Active principal application

#### SKIN

Acne

Face lifting

Wrinkle treatment

Scar

Couperose

Body spots

Face spots

Restructuring skin