

LIPOZERO EXCEL

AESTHETIC ULTRASOUND | 7 PROGRAMS

LIPOZERO EXCEL IS PART OF THE GLOBUS HOME CARE LINE, AND HAS BEEN STUDIED TO REALIZE AT HOME BEAUTY CENTRE PROFESSIONAL STYLE TREATMENTS. THE NEW PORTABLE DESIGN AND THE GRAPHIC MENU HAVE BEEN DESIGNED FOR THE HOME-DOMESTIC USE AND MAKES IT USERFRIENDLY EVEN FOR NON EXPERIENCED USERS. IT ALLOWS TO TREAT THE MOST COMMON IMPERFECTION DUE TO CELLULITE AND FAT ACCUMULATION.

5 BEAUTY | 2 CELLULITE

TECHNICAL FEATURES

Display

Visible area size 2,6"

Outputs

1

Head

1 head Ø mm 42 - IP68

Power

1,2W/cm² ± 20%

Emission Frequency

3MHz

Emission

Continuous and pulsed

Power Supply

- Mains
- Rechargeable batteries (not included)

Size


mm 100x160x35

Weight

gr 454

EQUIPMENT

- 1 Carrying case
- 1 Ultrasound LIPOZERO EXCEL
- 1 Head da 3MHz
- 1 Gel flask
- 1 Power supply unit
- 1 Operating manual



CE



MAINS SUPPLY



RECHARGEABLE BATTERIES

GLOBUS
ITALIAN EXCELLENCE

LIPOZERO EXCEL

AESTHETIC ULTRASOUND | 7 PROGRAMS

PROGRAMS LIST

BEAUTY

Stetch marks

Drainage

Massage

Modelling

Active principles application

CELLU-LYSE

Chock

Cellulite