

MAGNUM 3500 PRO

PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

THE ONLY GLOBUS MAGNETOTHERAPY DEVICE WITH FOUR CHANNELS THAT ALLOWS TO TREAT DIFFERENT BODY PARTS OR SUBJECTS AT THE SAME TIME. MAGNUM 3500 PRO IS EQUIPPED WITH A WIDE RANGE OF ACCESSORIES FOR ALL TREATMENT NEEDS.

TECHNICAL SPECIFICATIONS

Display

Size of the visible area 2.6"

Channels

4

Intensity can be emitted with flexible diffuser

- 200 Gauss on average for channel
- 800 Gauss on average on the whole (for 4 channels)

Frequency

5-200 Hz

Duty Cycle

5-50%

Free programming

Possibility of adding up to 10 custom programs

2+2 function

Possibility of running 2 different programs at the same time, using 2 channels each

Power supply

- Mains
- Rechargeable batteries (not included)

Size | Weight


170x220x60 mm | 820 gr



STANDARD EQUIPMENT

Code G5434

- 1 Suitcase
- 1 MAGNUM 3500 PRO magnetotherapy device
- 2 "Flexible" solenoids, 300x110 mm
- 4 "Pocket Pro" soft solenoids, 155x95 mm
- 1 "Soft 4" mat, 250x250 mm
- 2 Elastic bands with Velcro 800x80 mm
- 1 Magnum Test
- 1 Power supply unit
- 1 Operating manual



EASY PROGRAM



2+2 MODE



LAST 10



PROGRAMMABLE



MULTI USER



RENTAL OPTION



RUN TIME EDIT PARAMETERS



MAINS SUPPLY

MAGNUM 3500 PRO

PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

PROGRAM LIST

ARTHROSIS/ARTHRITIS

Knee arthrosis

Cervical arthrosis

Arthritis

Arthrosis

Hand arthrosis

Hip arthrosis

Shoulder arthrosis

Column arthrosis

PAIN

Algodystrophy

Back pain

Joint pain

INFLAMMATIONS

Anti-inflammatory

Tendinitis

Epicondylitis

Epitrochleytis

Carpal tunnel

Shoulder tendinitis

ORTHOPEDICS

Post ACL surgery

Cartilage lesions

Knee prosthesis

Hip prosthesis

Fusion

Osteotomy

BONE PATOLOGIES

Fractures

Scaphoid fracture

Wrist fracture

Shoulder fracture

Collarbone fracture

Forearm fracture

Pelvic fracture

Tibial fracture

Femoral fracture

Delayed fracture consolidation

Fracture with plaster

Osteobecrisis

Osteoporosis

Pseudoarthrosis

Bone edemas

TRAUMAS/EDEMAS/WOUNDS

Ankle sprain

Edemas - Contusion

Ulcers

NIGHT PROGRAMS

Night program 15 Hz

Night program 50 Hz

Night program 75 Hz

DEFAULT PROGRAMS WITH DIFFERENT FREQUENCY

Multi-frequency program

Program 5 Hz

Program 10 Hz

Program 15 Hz

Program 20 Hz

Program 25 Hz

Program 30 Hz

Program 35 Hz

Program 40 Hz

Program 45 Hz

Program 50 Hz

Program 55 Hz

Program 60 Hz

Program 65 Hz

Program 70 Hz

Program 75 Hz

Program 80 Hz

Program 85 Hz

Program 90 Hz

Program 95 Hz

Program 100 Hz

Program 120 Hz

Program 140 Hz

Program 160 Hz

Program 180 Hz

Program 200 Hz