

STIMVET 2000

4 CHANNELS ELECTROTHERAPY

STIMVET 2000 IS A PORTABLE ELECTROSTIMULATOR FOR THE MUSCULAR REHABILITATION AND FOR THE ANTALGIC TREATMENT OF CATS, DOGS AND HORSES. THE DEVICE IS ALSO SUITABLE FOR ATHLETE OR WORKING DOGS WHO HAVE TO INCREASE THEIR ATHLETIC SKILLS. STIMVET 2000 ALSO INCLUDES SOME MICROCURRENT PROGRAMS USEFUL FOR WOUND CICATRIZATION.


SPECIFIC PROGRAMS FOR DOGS, CATS AND HORSES

TECHNICAL FEATURES

Display
Backlit, 2.6"
Channels
4 independent (8 electrodes)
Frequency
1-150Hz
Pulse amplitude
50-450µs
Intensity
100mA per channel
Power supply
- Mains
- Battery pack
Size
100x160x35 mm
Weight
440 gr

STANDARD EQUIPMENT

- 1 Carrying case
- 1 STIMVET 2000 device
- 4 Elastic bands
- 2 Grey cables for microcurrents
- 4 Coloured cables for electrode connection
- 4 Silicone squared electrodes
- 4 Silicone rectangular electrodes
- 1 260 ml gel bottle
- 1 Power supply unit
- 1 Operating manual



- TENS** TENS
- NMES** NMES
- MENS** MENS
- DEN** DENERVATED
- RUN TIME** RUN TIME
- Last 10** LAST 10
- FAVORITES** FAVORITES FUNCTION
- MULTI USER** MULTI USER
- PRO** PROGRAMMABLE
- STIM LOCK** STIM LOCK
- 2+2** 2+2 MODE
- SyS** SYNCRO STIM
- BATTERY PACK** BATTERY PACK
- MAINS POWER SUPPLY** MAINS POWER SUPPLY

STIMVET 2000

4 CHANNELS ELECTROTHERAPY

Rev. En_03_2021

PROGRAM LIST

DOG

REHABILITATION

Recovery after ACL surgery

Contractures

Atrophy recovery 1 - 2

Muscle rebuilding 1 - 2

Strengthening 1 - 2

Motor point pen

Agonist / Antagonist

TENS

Endorphinic TENS

Conventional antalgic TENS

Anti-accommodation modulated TENS

Post-surgical pain

Fracture pain

Osteoarthritis

Wobbler syndrome

Discopathy

Spondylarthritis

Lumbosacral stenosis

SPORT

Forelimb maximum strength

Hind limbs maximum strength

Forelimb endurance strength

Hind limbs endurance strength

Forelimb explosive strength

Hind limb explosive strength

Decontracting

Cool down

DENERVATED

Triangular 1 - 2 - 3

Rectangular 1 - 2 - 3

Trapezoidal 1 - 2 - 3

MICROCURRENTS

Osteogenesis

Wound cicatrization

Cartilage

CAT

REHABILITATION

Recovery after ACL surgery

Contractures

Atrophy recovery 1 - 2

Muscle rebuilding 1 - 2

Strengthening 1 - 2

Motor point pen

Agonist / Antagonist

TENS

Endorphinic TENS

Conventional antalgic TENS

Anti-accommodation modulated TENS

Post-surgical pain

Fracture pain

Osteoarthritis

Discopathy

Spondylarthritis

Lumbosacral stenosis

DENERVATED

Triangular 1 - 2 - 3

Rectangular 1 - 2 - 3

Trapezoidal 1 - 2 - 3

MICROCURRENTS

Osteogenesis

Wound cicatrization

Cartilage

HORSE

REHABILITATION

Recovery after ACL surgery

Contractures

Atrophy recovery 1

Motor point pen

Agonist / Antagonist

TENS

Endorphinic TENS

Conventional antalgic TENS

Post-surgical pain

Fracture pain

Osteoarthritis

Discopathy

Spondylarthritis

MICROCURRENTS

Osteogenesis

Wound cicatrization

Cartilage