

STIMVET 4000

4 CHANNELS ELECTROTHERAPY

STIMVET 4000 IS A PROFESSIONAL ELECTROTHERAPY DEVICE FOR VETERINARY REHABILITATION. THE DEVICE INCLUDES CURRENTS SUITABLE FOR MUSCLE RECOVERY IN CASE OF ATROPHY AND TENS CURRENTS FOR THE ANALGIC TREATMENT OF CHRONIC AND ACUTE PAIN. IN ADDITION, THERE ARE CURRENTS TO TREAT DENERVATED MUSCLES AND MICROCURRENTS, MAINLY USED FOR CARTILAGE LESIONS AND WOUNDS. THE DEVICE IS FULLY EQUIPPED WITH DEFAULT PROGRAMS FOR THE TREATMENT OF CATS, DOGS AND HORSES.

SPECIFIC PROGRAMS FOR DOGS, CATS AND HORSES

TECHNICAL FEATURES

Display

Backlit, 2.6"

Channels

4 independent (8 electrodes)

Frequency

1-150Hz

Pulse amplitude

50-450µs

Intensity

100mA per channel

Power supply

- Mains

- Battery pack

Size



170x220x60 mm

Weight

1370 gr

STANDARD EQUIPMENT

- 1 Carrying case
- 1 STIMVET 4000 device
- 4 Elastic bands
- 2 Grey cables for microcurrents and ionophoresis
- 4 Coloured cables for electrode connection
- 4 Silicone squared electrodes
- 4 Silicone rectangular electrodes
- 4 Sponges for ionophoresis
- 1 260 ml gel bottle
- 1 Power supply unit
- 1 Operating manual



STIMVET 4000

4 CHANNELS ELECTROTHERAPY

Rev. En_03_2021

PROGRAM LIST

DOG

REHABILITATION

Recovery after ACL surgery

Contractures

Atrophy recovery 1 - 2

Muscle rebuilding 1 - 2

Strengthening 1 - 2

Motor point pen

Agonist / Antagonist

TENS

Endorphinic TENS

Conventional antalgic TENS

Anti-accommodation modulated TENS

Post-surgical pain

Fracture pain

Osteoarthritis

Wobbler syndrome

Discopathy

Spondylarthritis

Lumbosacral stenosis

SPORT

Forelimb maximum strength

Hind limbs maximum strength

Forelimb endurance strength

Hind limbs endurance strength

Forelimb explosive strength

Hind limb explosive strength

Decontracting

Cool down

DENERVATED

Triangular 1 - 2 - 3

Rectangular 1 - 2 - 3

Trapezoidal 1 - 2 - 3

MICROCURRENTS

Osteogenesis

Wound cicatrization

Cartilage

IONOPHORESIS

Ionophoresis

CAT

REHABILITATION

Recovery after ACL surgery

Contractures

Atrophy recovery 1 - 2

Muscle rebuilding 1 - 2

Strengthening 1 - 2

Motor point pen

Agonist / Antagonist

TENS

Endorphinic TENS

Conventional antalgic TENS

Anti-accommodation modulated TENS

Post-surgical pain

Fracture pain

Osteoarthritis

Discopathy

Spondylarthritis

Lumbosacral stenosis

DENERVATED

Triangular 1 - 2 - 3

Rectangular 1 - 2 - 3

Trapezoidal 1 - 2 - 3

MICROCURRENTS

Osteogenesis

Wound cicatrization

Cartilage

IONOPHORESIS

Ionophoresis

HORSE

REHABILITATION

Recovery after ACL surgery

Contractures

Atrophy recovery 1

Motor point pen

Agonist / Antagonist

TENS

Endorphinic TENS

Conventional antalgic TENS

Post-surgical pain

Fracture pain

Osteoarthritis

Discopathy

Spondylarthritis

MICROCURRENTS

Osteogenesis

Wound cicatrization

Cartilage

IONOPHORESIS

Ionophoresis