# **GENESY 1500**

# ELECTROTHERAPY 4 CHANNELS | 414 PROGRAMS

PROFESSIONAL FOUR-CHANNEL PRODUCT USED IN THE MEDICAL FIELD, COMPLETE OF ALL THE CURRENTS FOR PHYSIOTERAPY AND REHABILITATION THERAPY.

92 ELECTROTHERAPY MENÙ (EMS/RUSSIAN) | 84 ACTION NOW | 54 SERIAL SEQUENTIAL STIMULATION | 53 SPORT 33 PAIN | 29 FITNESS | 23 MICROCURRENTS | 16 BEAUTY | 11 REHAB | 9 DENERVATED | 6 INTERFERENTIALS 3 INCONTINENCE | 1 IONOPHORESIS

#### **TECHNICAL FEATURES**

### Display

Visible area size 2,6"

#### Channels

4 indipendent (8 electrodes)

#### Frequency

0,3-150 Hz \*

#### Pulse amplitude

40-450 us \*

#### **Power**

0-120 mA per channel

#### Power supply

- Mains
- Rechargeable Batteries

#### Size

100x160x35 mm

#### Weight

450 gr

\* According to the type of current

## **EQUIPMENT**

- · 1 Bag
- · 1 GENESY 1500 stimulator
- · 4 Cables for electrode connection
- · 2 Cables for microcurrents and ionophoresis
- · 4 Self-adhesive square electrodes
- · 4 Self-adhesive rectangular electrodes
- · 1 Power supply unit
- · 1 Operating manual


























35







**MULTI USER** 

**FAVORITES** 





















# **GENESY 1500**

# ELECTROTHERAPY 4 CHANNELS | 414 PROGRAMS

# **PROGRAMS LIST**

REHAB	Modulated antalgic tens	SPORT
Swollen ankles	IONOPHORESIS	Capillarization
Atrophy recovery	lonophoresis	Warm-up
Hemiplegy-upper limbs	INCONTINENCE	Pre-competition w
Hemiplegy-lower limbs	Mixed incontinence	Active recovery
Recovery after ACL surgery	Stress incontinence	Maximum strength
Functional recovery	URGE incontinence	Endurance strengt
Ankle re-education	MICROCURRENTS	Explosive strength
Leg re-education	Epicondylitis	Aerobic enduranc
Muscle Spasms	Scapulohumeral periarthritis	Reactivity
Shoulder subluxation prevention	Muscle restoration	Post-competition r
Vastus medialis reinforcement	Contusion	Decontracting
Motor point pen	Edema	Hypertrophy
Agonist/Antagonist	Skin ulcera	FITNESS
Muscle reinforcement	Sciatica	Firming
PAIN	Lumbago	Bio-Pulse firming
High frequency antalgic tens	Brachial neuralgia	Sculpting
Low frequency antalgic tens	Acute pain	Bio-Pulse sculpting
Endorphinic tens	Articular pain	Toning
Knee pain	Stiff neck	Mass building
Post-surgery pain	Whiplash	Body sculpting
Chronic pain	Cervical spondylosis	 Definition
Shoulder pain	Shoulder sprain	Jogging
Carpal tunnel	Carpal tunnel	Anaerobic fitness
Trapezius pain	Knee sprain	BEAUTY
Knee osteoarthritis	Knee osteoarthritis	 Drainage
Rotator cuff tendinitis	Ankle sprain	 Lipolysis
Muscle pain	Achille tendon inflammation	Toning massage
Nerve compression	Patella tendon inflammation	Skin tone improver
Spinal compression	Rotator cuff inflammation	Connective masso
Osteoarthritis	Tendon inflammation	Post-pregnancy di
Spinal osteoarthritis	DENERVATED	Post-pregnancy lip
Spinal osteoporosis	Triangular 1	Post-pregnancy fir
Hip osteoarthritis	Triangular 2	
Ankle osteoarthritis	Triangular 3	The "3S" program
Muscle tendon injury	Trapezoidal 1	combinations.  The "3S" programs of
Post-surgical pain	Trapezoidal 2	activation of the char the channels 1 and 2.
Sciatica	Trapezoidal 3	The Serial Sequent stimulate the muscul
Chronic lumbago	Rectangular 1	to the differentiate muscular groups invo
Cervical pain	Rectangular 2	
Acute pain post inguinal hernia	Rectangular 3	ACTION NOW The Action Now progr
Bursitis-tendinitis	INTERFERENTIALS	combinations.
Muscle injuries	Knee osteoarthritis	ELECTROTHERAPY
Bone fractures	Acute lumbago	The Electrotherapy i EMS and 44 Russian (
Whiplash	Cervical pain	
Menstrual pain	Chronic lumbago	

#### **SPORT**

Capillarization Warm-up Pre-competition warm-up Active recovery Maximum strength Endurance strength Explosive strength Aerobic endurance Reactivity Post-competition recovery Decontracting Hypertrophy **FITNESS** Firming Bio-Pulse firming Sculpting Bio-Pulse sculpting Toning Mass building Body sculpting Definition

## **BEAUTY**

Drainaae Lipolysis Toning massage Skin tone improvement Connective massage Post-pregnancy drainage Post-pregnancy lipolysis Post-pregnancy firming

# **ERIAL SEQUENTIAL STIMULATION**

The "3S" program list includes 54 parameter combinations.

The "3S" programs are characterized by a delay activation of the channels 3 and 4 compared with

the channels 1 and 2.
The Serial Sequential Stimulation permits to stimulate the musculature in kinetic chain thanks to the differentiated activation times of the muscular groups involved.

### CTION NOW

The Action Now program list includes 84 parameter combinations.


#### LECTROTHERAPY

The Electrotherapy includes 48 combinations of EMS and 44 Russian (KOTZ) parameters.



Conventional antalgic tens

**Epicondylitis** 



Frozen shoulder

Post-surgical pain



