

# GENESY 3000

## ELECTROTHERAPY 4 CHANNELS | 423 PROGRAMS

TOP PROFESSIONAL FOUR-CHANNEL PRODUCT, USED IN THE MEDICAL FIELD FOR THE TREATMENT OF THE MOST COMMON PATHOLOGIES OF THE MUSCLE-SKELETICAL SYSTEM.

236 ELECTROTHERAPY MENU (EMS/RUSSIAN) | 53 SPORT | 33 PAIN | 29 FITNESS | 23 MICROCURRENTS | 16 BEAUTY | 11 REHAB | 9 DENERVATED | 9 INTERFERENTIALS | 3 INCONTINENCE | 1 IONOPHORESIS


### TECHNICAL FEATURES

<b>Display</b>
Visible area size 2,6"
<b>Channels</b>
4 independent (8 electrodes)
<b>Frequency</b>
0,3-150 Hz *
<b>Pulse amplitude</b>
40-450 µs *
<b>Power</b>
0-120 mA per channel
<b>Power supply</b>
- Mains - Rechargeable batteries
<b>Size</b>
170x220x60 mm
<b>Weight</b>
1,120 Kg

\* According to the type of current

### EQUIPMENT

- 1 Bag
- 1 GENESY 3000 stimulator
- 4 Cables for electrode connection
- 2 Cables for microcurrents and ionophoresis
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 2 Iono electrodes + 2 elastic bands
- 1 Power supply unit
- 1 Operating manual



**CE**  
0476



- |              |                        |
|--------------|------------------------|
| TENS         | STIM LOCK              |
| EMS          | FUNCTION 2+2           |
| IONO         | MY TRAINER             |
| RUN TIME     | AUTO STIM              |
| LAST 10      | SYNCRO STIM            |
| FAVORITES    | WORK TIME              |
| MULTI USER   | MAIN SUPPLY            |
| PROGRAMMABLE | RECHARGEABLE BATTERIES |

**GLOBUS**  
ITALIAN EXCELLENCE

# GENESY 3000

## ELECTROTHERAPY 4 CHANNELS | 423 PROGRAMS

### PROGRAMS LIST

Rev. En\_05\_2022

#### REHAB

Swollen ankles

Atrophy recovery

Hemiplegy-upper limbs

Hemiplegy-lower limbs

Recovery after ACL surgery

Functional recovery

Ankle re-education

Leg re-education

Muscle Spasms

Shoulder subluxation prevention

Vastus medialis reinforcement

Motor point pen

Agonist/Antagonist

Muscle reinforcement

#### PAIN

High-frequency antalgic tens

Low frequency antalgic tens

Endorphinic tens

Knee pain

Post-surgical pain

Chronic pain

Shoulder pain

Carpal tunnel

Trapezius pain

Osteoarthritis

Knee osteoarthritis

Rotator cuff tendinitis

Muscle pain

Nerve compression

Spinal compression pain

Spinal osteoarthritis

Spinal osteoporosis

Hip osteoarthritis

Ankle osteoarthritis

Muscle tendon injury pain

Post-surgical pain

Sciatica

Chronic lumbago

Cervical pain

Acute pain inguinal hernia

Bursitis-tendinitis

Muscle injuries

Bone fractures

Whiplash

Epicondylitis

Menstrual pain

Conventional antalgic tens

Modulated antalgic tens

#### IONOPHORESIS

Ionophoresis

#### INCONTINENCE

Mixed incontinence

Stress incontinence

URGE incontinence

#### MICROCURRENTS

Epicondylitis

Scapulothoracic periarthritis

Muscle restoration

Contusion

Edema

Skin ulcer

Sciatica

Lumbago

Brachial neuralgia

Acute pain

Articular pain

Stiff neck

Whiplash

Cervical spondylosis

Shoulder sprain

Carpal tunnel

Knee sprain

Knee osteoarthritis

Ankle sprain

Achilles tendon inflammation

Patella tendon inflammation

Rotator cuff inflammation

Tendon inflammation

#### DENERVATED

Triangular 1

Triangular 2

Triangular 3

Trapezoidal 1

Trapezoidal 2

Trapezoidal 3

Rectangular 1

Rectangular 2

Rectangular 3

#### INTERFERENTIALS

Knee osteoarthritis

Acute lumbago

Cervical pain

Chronic lumbago

Frozen shoulder

Knee post-surgery pain

Pain vertebral fractures

Knee post-surgery edema

Fibromyalgia

#### SPORT

Capillarization

Warm-up

Pre-competition warm-up

Active recovery

Maximum strength

Endurance strength

Explosive strength

Aerobic endurance

Reactivity

Post-competition recovery

Decontracting

Hypertrophy

#### FITNESS

Firming

Bio-Pulse firming

Sculpting

Bio-Pulse sculpting

Toning

Mass building

Body sculpting

Definition

Jogging

Anaerobic fitness

Aerobic fitness

Cramp prevention

#### BEAUTY

Drainage

Lipolysis

Toning massage

Skin tone improvement

Connective massage

Post-pregnancy drainage

Post-pregnancy lipolysis

Post-pregnancy firming

#### ELECTROTHERAPY

The Electrotherapy program list includes 192 EMS and 44 Russian(KOTZ) parameters combinations