

MAGNETOVET 200

MAGNETOTHERAPY

MAGNETOVET 200 IS A HIGH INTENSITY PULSED ELECTROMAGNETIC FIELD MAGNETOTHERAPY EQUIPMENT. THE IMPULSES PRODUCED HAVE A THERAPEUTIC AND REGENERATIVE EFFECT ON THE MAIN OSTEOARTICULAR PATHOLOGIES THAT AFFECT PETS. FOR ITS SIMPLICITY OF USE MAGNETOVET 200 IS PARTICULARLY SUITABLE ALSO FOR HOME TREATMENT OF CHRONIC PAIN.

SPECIFIC PROGRAMS FOR DOGS, CATS AND HORSES

TECHNICAL FEATURES

Display

Backlit, 2.6"

Channels

2

Emitted intensity

200 average Gauss per channel with flex solenoids

Frequency

Adjustable in programming from 5 to 200Hz

Duty Cycle

From 5% to 50% in steps of 5

Power supply

- Mains
- Battery pack

Size



100x160x35 mm

Weight

280 gr

STANDARD EQUIPMENT

- 1 Carrying case
- 1 MAGNETOVET 200 device
- 1 Flexible solenoids 300x100 mm
- 1 Power supply unit
- 1 Operating manual



MAGNETOVET 200

MAGNETOTHERAPY

Rev. En_05_2024

PROGRAM LIST

DOG / CAT

Fractures

Pseudoarthrosis

Wounds / Ulcers

Hip osteoarthritis

Knee osteoarthritis

Elbow osteoarthritis

Shoulder osteoarthritis

Spondylarthrosis

Luxation

Cicatrization

Edema

Peripheral paralyses

Inflammation pain

HORSE

Osteoarthritis

Wounds / Ulcers

Discopathy

Edema