

RF CLINIC BODY

RADIO FREQUENCY | 10 PROGRAMS

RF CLINIC BODY IS AN INNOVATIVE PROFESSIONAL DEVICE PRODUCING RADIO FREQUENCY WAVES FOR BEAUTY APPLICATIONS AGAINST THE AGING OF FACE AND BODY SKIN. RF CLINIC BODY ENHANCES THE VASCULARIZATION AND THE REGENERATION OF THE SKIN TISSUES, ELIMINATING THE IMPERFECTIONS DUE TO THE SUN EXPOSURE AND TO THE SKIN AGING.

6 SKIN | 4 FACE

TECHNICAL FEATURES

Display

Visible area size 2,6"

Outputs

1

Head

1 head RF Clinic 100

Power

200W (50Ω)

Emission frequency

From 470KHz to 2MHz

Treatment methods

Resistive

Power supply

- Mains
- Rechargeable batteries (not included)

Size


170x220x60 mm

Weight

920 gr

EQUIPMENT

- 1 Carrying case
- 1 Radio frequency unit RF CLINIC BODY
- 1 Head RF Clinic 100
- 1 Cream can 250 ml
- 1 Power supply unit
- 1 Operating manual



CE



MAINS SUPPLY



RECHARGEABLE BATTERIES

GLOBUS
ITALIAN EXCELLENCE

RF CLINIC BODY

RADIO FREQUENCY | 10 PROGRAMS

PROGRAMS LIST

Rev. EN_03_2024

FACE

- Face wrinkles
- Eye wrinkles
- Face acne
- Face microlifting

SKIN

- Body wrinkles
- Sensitive skin wrinkles
- Cleavage micro-lifting
- Stretch marks
- Skin elasticity
- Bioskin collagen