

RF CLINIC PRO

RADIO FREQUENCY | 30 PROGRAMS

RF CLINIC PRO IS AN INNOVATIVE PROFESSIONAL DEVICE CREATING RADIO FREQUENCY WAVES FOR BEAUTY APPLICATIONS AGAINST THE AGING OF FACE AND BODY SKIN. THANKS TO THE BIPOLAR HANDPIECE OF 50MM IT IS POSSIBLE TO MAKE FIRING PROGRAMS AND FOR THE CELLULITE TREATMENT.

13 SKIN | 10 CELLULITE | 7 FACE

TECHNICAL FEATURES

Display

Visible area size 2,6"

Outputs

1

Head

- 1 head Ø 50 mm
- 1 head RF Clinic 100

Power

200W (50Ω)

Emission frequency

from 470KHz to 2MHz

Treatment methods

Resistive

Power supply

- Mains
- Rechargeable batteries (not included)

Size



170x220x60 mm

Weight

920 gr

EQUIPMENT

- 1 Carrying case
- 1 Radio frequency unit RF CLINIC PRO
- 1 Head RF Clinic 500
- 1 Head RF Clinic 100
- 1 Cream can 250 ml
- 1 Power supply unit
- 1 Operating manual



CE



MAINS SUPPLY



RECHARGEABLE BATTERIES

GLOBUS
ITALIAN EXCELLENCE

RF CLINIC PRO

RADIO FREQUENCY | 30 PROGRAMS

PROGRAMS LIST

FACE

Face wrinkles

Eye wrinkles

Face acne

Face micro-ifting

Couperosis

Face pimples

Anti aging

SKIN

Skin acne

Body wrinkles

Sensitive skin wrinkles

Cleavage microlifting

Stretch marks

Skin elasticity

Bioskin collagen

Scars

Body pimples

Skin restructure

Skin smoothing

Firming

Massage

CELLULITE

Lipolysis

Cellulitis

Edematous cellulite

Fibrous cellulite

Flabby thighs

Orange skin cellulite

Man abdomen adiposity

Man hip adiposity

Lower limb man abdomen adiposity

Upper limb man abdomen adiposity