LIPOZERO G800

AESTHETIC ULTRASOUND | 30 PROGRAMS

DEVICE FOR PROFESSIONAL AESTHETICS, SUITABLE FOR FIGHTING THE IMPERFECTIONS OF CELLULITE AND LOCALIZED ADIPOSITY. AND TO IMPROVE SKIN APPEARANCE. THE POSSIBILITY OF PROVIDING DIFFERENT FREQUENCIES WILL GUARANTEE A SPECIFIC AND EFFECTIVE TREATMENT FOR EVERY GOAL.

10 MF-PLUS | 10 BEAUTY | 10 SKIN

TECHNICAL SPECIFICATIONS

Display

Size of the visible area 2.6"

Outputs

1

Handpiece

1 handpiece Ø mm 42 - IP68

Power

- 3 W/cm² ± 20% a 1.3 MHz e 3 MHz

- 1,5 W/cm² ± 20% a 1 MHz e 800 KHz

Emitted frequency

3 MHz - 1,3 MHz - 1 MHz - 800 KHz

Issue

Continuous and pulsed

Power supply

- Mains

- Rechargeable battery (not included)

Size

170x220x60 mm


Weight

920 gr



STANDARD EQUIPMENT

- · 1 Suitcase
- ·1 Ultrasound LIPOZERO G800
- ·1 Multifrequency handpiece, Ø 42 mm IP68
- · 1 Gel bottle
- ·1 Power supply
- · 1 Operating manual





CE

In compliance with the Interministerial Decree n. 206/2015, the device is suitable for use in beauty centers.



MAINS SUPPLY

RECHARGEABLE BATTERIES



LIPOZERO G800

AESTHETIC ULTRASOUND | 30 PROGRAMS

PROGRAM LIST

MF-PLUS

MF-PLUS
Lipolysis
Cellulite
Edematous cellulite
Fibrous cellulite
Flabby thighs
Orange peel cellulite
Male abdomen adiposity
Male hip adiposity
Upper limb male adiposity
Lower limb male adiposity
BEAUTY
Firming
Massage
Sculpting
Abdomen sculpting
Leg sculptin
Drainage
Arm drainage
Leg drainage
Stretch marks
Active principle application
SKIN
Acne
Face lifting
Wrinkle treatment
Scars
Rosacea
Face pimples
Body pimples
Skin restructuring

Skin smoothing

Photoaging





