LIPOZERO EXCEL

AESTHETIC ULTRASOUND | 7 PROGRAMS

LIPOZERO EXCEL IS PART OF THE GLOBUS HOME CARE LINE, AND HAS BEEN STUDIED TO REALIZE AT HOME BEAUTY CENTRE PROFESSIONAL STYLE TREATMENTS. THE NEW PORTABLE DESIGN AND THE GRAPHIC MENU HAVE BEEN DESIGNED FOR THE HOME-DOMESTIC USE AND MAKES IT USERFRIENDLY EVEN FOR NON EXPERIENCED USERS. IT ALLOWS TO TREAT THE MOST COMMON IMPERFECTION DUE TO CELLULITE AND FAT ACCUMULATION.

5 BEAUTY | 2 CELLULITE

TECHNICAL FEATURES

Display

Visible area size 2,6"

Outputs

Head

1 head Ø mm 42 - IP68

Power

1,2 W/cm² ± 20%

Emission Frequency

3 MHz

Emission

Continuous and pulsed

Power Supply

- Mains
- Rechargeable batteries (not included)

100x160x35 mm

Weight

454 gr

EQUIPMENT

- · 1 Carrying case
- · 1 Ultrasound LIPOZERO EXCEL
- · 1 3 MHz handpiece
- · 1 Gel flask
- · 1 Power supply unit
- · 1 Operating manual



















ALIMENTAZIONE A RETE



BATTERIA RICARICABILE



Rev. En_02_2022

LIPOZERO EXCEL

AESTHETIC ULTRASOUND | 7 PROGRAMS

PROGRAMS LIST

Stetch marks

Drainage

Massage

Modelling

Active principles application

CELLU-LYSE

Chock

Cellulite





