

LIPOZERO G39

AESTHETIC ULTRASOUND | 49 PROGRAMS

LIPOZERO G39 IS THE PROFESSIONAL LIPOLYTIC ULTRASOUND THAT FIGHTS CELLULITE AND ELIMINATES EXCESS FAT THANKS TO PROGRAMS CONCEIVED TO INTERVENE AT DIFFERENT FREQUENCIES AND ENSURE MAXIMUM EFFECTIVENESS IN EVERY BEAUTY TREATMENT.

10 MF-PLUS | 10 BEAUTY | 10 SKIN | 6 CAVITATION CAVIMAX 150 kHz

5 ADIPOSITY MOBILIZATION NorAd-526 | 4 MALE ADIPOSITY | 4 ULTRACAV 39 kHz

TECHNICAL SPECIFICATIONS

Display

Size of the visible area 2.6"

Outputs

1

Handpiece

- 1 handpiece Ø mm 50 - IP68

- 1 handpiece Ø mm 42 - IP68

Power

3 W/cm² ± 20%

Emitted frequency

3 MHz - 1 MHz - 526 kHz - 150 kHz - 39 kHz

Emission

Continuous and pulsed

Power supply

- Mains

- Rechargeable battery (not included)

Size

170x220x60 mm

Weight

920 gr

STANDARD EQUIPMENT

· 1 Suitcase

· 1 Ultrasound LIPOZERO G39


· 1 Multifrequency handpiece

· 1 Handpiece to use with 39 kHz

· 1 Gel bottle

· 1 Power supply

· 1 Operating manual



MAINS SUPPLY



RECHARGEABLE BATTERIES

LIPOZERO G39

AESTHETIC ULTRASOUND | 49 PROGRAMS

PROGRAM LIST

ULTRACAV 39 KHZ

Total body 3-6 cm

Total body >6 cm

Thighs/gluteus 3-6 cm

Thighs/gluteus >6 cm

CAVIMAX 150 KHZ

Lipolysis

Abdominal area lipolysis

Cellulite

Edematous cellulite

Fibrous cellulite

Orange peel cellulite

MF-PLUS

Lipolysis

Cellulite

Edematous cellulite

Fibrous cellulite

Flabby thighs

Orange peel cellulite

Male abdomen adiposity

Male hip adiposity

Upper limb male adiposity

Lower limb male adiposity

ADIPOSIY MOBILIZATION NorAd-526

Lipolysis NorAd-526

Cellulite NorAd-526

Adiposity NorAd-526

NorAd-526 1

NorAd-526 2

MALE ADIPOSIY

Male abdomen adiposity

Male hip adiposity

Upper limb male adip

Lower limb male adiposity

BEAUTY

Firming

Massage

Sculpting

Abdomen sculpting

Leg sculptin

Drainage

Arm drainage

Leg drainage

Stretch marks

Active principle application

SKIN

Acne

Face lifting

Wrinkle treatment

Scars

Rosacea

Face pimples

Body pimples

Skin restructuring

Skin smoothing

Photoaging