

PODCARE 12.0 PRO

LASER THERAPY


TOP OF THE RANGE AMONG THE LASER DEVICES FOR PODOLOGY. WITH ITS 12W OF MAXIMUM POWER, IT ALLOWS TO PERFORM DIFFERENT PODOLOGICAL AND MUSCULOSKELETAL TREATMENTS OPTIMIZING TIME AND GUARANTEEING A DEEP ACTION.

TECHNICAL SPECIFICATIONS

Display
2.6", backlit
Laser source
Class 4 GaAlAs diode
Maximum power
12W \pm 20%
Compatible with a wavelength of
808nm
Issue
Pulsed and continuous
Target area dimensions
11 mm optic 22 mm opti
Frequency
Up to 10,000Hz
Power supply
Mains and rechargeable batteries
Size and weight
220x170x60 mm - 1.12 Kg

STANDARD EQUIPMENT

- 1 Suitcase
- 1 PODCARE 12.0 PRO device
- 1 808nm laser handpiece with 11mm and 22mm optics
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTIONS

- Cod. G5995 HP 980nm handpiece
- Cod. G5996 1064nm handpiece, maximum power 3W



CE
0476



FAVORITES



PROGRAMMABLE



BATTERY PACK



MAINS POWER SUPPLY

GLOBUS
ITALIAN EXCELLENCE

PODCARE 12.0 PRO

LASER THERAPY

MEDICAL PROGRAM LIST

808nm WAVELENGTH PROGRAMS

MUSCULOSKELETAL

GENERAL

Arthrosis

Rheumatoid arthritis

Muscle injuries

Muscle contracture

Radiculopathies

Tendinopathies

Neuralgia

Posttherpetic neuralgia

Sudeck disease

Bursitis

Edema

HEAD - NECK

Arthrosis

Radiculopathies

Myofascial pain

Chronic cervical

Temporo-mandibular pain

Rheumatoid arthritis

Trigeminal neuralgia

Cervical herniation

ARM - HAND

Arthrosis

Rheumatoid arthritis

Raynaud syndrome

Finger pain

Epicondylitis

Carpal tunnel

Wrist pain

Tendinitis

KNEE

Arthrosis

Rheumatoid arthritis

Meniscus pain

Patello-femoral

Ligament injury

Gonalgia

Sprain

Chondropathy

Patello-femoral tendinitis

FOOT - ANKLE

Arthrosis

Rheumatoid arthritis

Onychocryptosis

Onychomycosis

Metatarsalgia

Ankle pain

Plantar fasciitis

Extensor tendinitis

Distal tendinopathies

Achilles tendon

Tibial tendinitis

Ankle tendinitis

Peroneal tendonitis

Sprain

Edema

Morton's neuroma

Heel pain

Raynaud syndrome

BACK

Vertebral arthrosis

Acute pain

Lumbago

Radiculopathies

Rheumatoid arthritis

SHOULDER

Myofascial pain

Arthrosis

Adhesive capsulitis

Bursitis

Subacromial impingement

Tendinitis

Luxations

Tendon injuries (cuff)

Rheumatoid arthritis

Edema

DERMATOLOGY

Ulcers

Cicatrization_ulcers

Wounds

980nm WAVELENGTH PROGRAMS

Antalgic

Cervical Pain

Lumbar pain

Bursitis

Baker's cyst

Cruralgia

De Quervain syndrome

Sprains

Contusions

Acute pain

Chronic pain

Edemas

Hematomas

Extended hematomas

Epicondylitis

Epitrochleitis

Discal hernia

Fibromyalgia

Inflammations

Ligament injury

Muscle injuries

Meniscus pain

Synovitis

Muscle strains

Tendinopathy

Carpal Tunnel