

PODCARE 6.0

LASER THERAPY

PODCARE 6.0 IS A LIGHTWEIGHT AND PORTABLE DEVICE, SUITABLE TO ALWAYS BE CARRIED ALONG TO PERFORM TREATMENTS IN MULTIPLE CLINICS. IT IS EQUIPPED WITH A 808NM DIODE WITH A MAXIMUM POWER OF 6W. THE DEVICE IS IDEAL FOR PODOLOGICAL APPLICATIONS ON THE FOOT AND THE ANKLE, BUT IT ALSO HAS A SERIES OF SPECIFIC PROGRAMS FOR ULCERS, WOUNDS AND MUSCULOSKELETAL DISORDERS.

TECHNICAL SPECIFICATIONS

Display

2.6", backlit

Laser source

Class 4 GaAlAs diode

Maximum power

6W \pm 20%

Compatible with a wavelength of

808nm

Issue

Pulsed and continuous

Target area dimensions

11 mm optic

22 mm optic (optional)

Frequency

Up to 10,000Hz

Power supply


Mains and rechargeable batteries

Size and weight

100x160x35 mm - 450 g

STANDARD EQUIPMENT

- 1 Suitcase
- 1 PODCARE 6.0 device
- 1 808nm laser handpiece with 11mm optic
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTIONS

- Cod. G5994 MMP 980nm handpiece
- Cod. G5996 1064nm handpiece, maximum power 3W
- Cod. G6092 22 mm optic



CE
0476



FAVORITES



PROGRAMMABLE



BATTERY PACK



MAINS POWER SUPPLY

GLOBUS
ITALIAN EXCELLENCE

PODCARE 6.0

LASER THERAPY

MEDICAL PROGRAM LIST

808nm WAVELENGTH PROGRAMS

MUSCULOSKELETAL

GENERAL

Arthrosis

Rheumatoid arthritis

Muscle injuries

Muscle contracture

Fibromyalgia

Radiculopathies

Tendinopathies

Neuralgia

Bursitis

Edema

HEAD - NECK

Arthrosis

Radiculopathies

Chronic cervical

Temporo-mandibular pain

Rheumatoid arthritis

ARM - HAND

Arthrosis

Rheumatoid arthritis

Raynaud syndrome

Finger pain

Epicondylitis

Carpal tunnel

Wrist pain

Tendinitis

KNEE

Arthrosis

Rheumatoid arthritis

Meniscus pain

Bursitis

Patello-femoral

Ligament injury

Sprain

Chondropathy

Patello-femoral tendinitis

FOOT - ANKLE

Arthrosis

Rheumatoid arthritis

Onychocryptosis

Onychomycosis

Metatarsalgia

Ankle pain

Plantar fasciitis

Extensor tendinitis

Distal tendinopathies

Achilles tendon

Ankle tendinitis

Sprain

Morton's neuroma

Heel pain

Raynaud syndrome

BACK

Vertebral arthrosis

Lumbago

Myofascial pain

Radiculopathies

Rheumatoid arthritis

SHOULDER

Myofascial pain

Arthrosis

Adhesive capsulitis

Bursitis

Subacromial impingement

Tendinitis

Luxations

Tendon injuries (cuff)

Rheumatoid arthritis

DERMATOLOGY

Ulcers

Cicatrization_ulcers

Wounds

980nm WAVELENGTH PROGRAMS

Antalgic

Cervical Pain

Lumbar pain

Bursitis

Baker's cyst

Cruralgia

De Quervain syndrome

Sprains

Contusions

Acute pain

Chronic pain

Edemas

Hematomas

Extended hematomas

Epicondylitis

Epitrochleitis

Discal hernia

Fibromyalgia

Inflammations

Ligament injury

Muscle injuries

Meniscus pain

Synovitis

Muscle strains

Tendinopathy

Carpal Tunnel