

# TECAR BEAUTY ELITE

## AESTHETIC TECARTHERAPY

TECAR BEAUTY IS A NON-INVASIVE METHOD THAT FIGHTS AGING AND SKIN LAXITY. THE ELECTROMAGNETIC WAVES PRODUCED BY THE DEVICE PENETRATE INTO THE DERMIS AND RELEASE HEAT IN THE TISSUES. THIS CAUSES THE SYNTHESIS OF NEW COLLAGEN AND ELASTIN, VISIBLY IMPROVING THE COMPACTNESS OF THE SKIN. IN ADDITION, THE HEAT FACILITATES MICROCIRCULATION, PROMOTING A BETTER OXYGENATION OF THE TISSUES. PROFESSIONAL DEVICE SUITABLE FOR AESTHETIC CENTERS.

### TECHNICAL SPECIFICATIONS

#### Display

Size 7" color touch screen

#### Power

50 W  $\pm$ 10% Capacitive

25 W  $\pm$ 10% Resistive

#### Emission frequencies

400 kHz, 448 kHz, 470 kHz, 700 kHz,

1000 kHz, 1200 kHz,  $\pm$  10%

#### Treatment methods

Capacitive and Resistive

#### Capacitive Electrodes

3  $\varnothing$  (30, 50, 70 mm)

#### Resistive Electrodes

3  $\varnothing$  (30, 50, 70 mm)

#### Mains power supply

100-240 V AC, 50-60 Hz, 1.8 A

#### Size

515x340x170 mm

#### Weight

11,5 Kg

### STANDARD EQUIPMENT

· Cod. G6243

· 1 TECAR BEAUTY ELITE device

· 1 Counter electrode (flexible plate)

· 1 Monopolar Resistive handpiece with 3 electrodes  
 $\varnothing$  (30, 50, 70 mm)

· 1 Monopolar Capacitive handpiece with 3 electrodes  
 $\varnothing$  (30, 50, 70 mm)

· 1 Bipolar handpiece  $\varnothing$  (30 mm)

· 1 Hands-free treatment kit

· 1 500 ml hyaluronic acid conductive face cream

· 1 1000 ml conductive cream for lymphatic-drainage  
with eascin, fucus and caffeine

· 1 Power supply cable

· 1 Operating manual













In compliance with Interministerial Decree no. 206/2015, the device is suitable for beauty centers.



# TECAR BEAUTY ELITE

## AESTHETIC TECARTHERAPY

### ESTHETIC ATLAS PROGRAM LIST

<b>FACE</b>	<b>GLUTEUS</b>	 TREATMENT ATLAS
Wrinkles	Stretch marks	 RESISTIVE
Eye contour	Firming	 CAPACITIVE
Acne scars	Scars	 MULTI-FREQUENCY
Lifting	Lymphatic drainage	 MYOFASCIALS
Rosacea	Fat reduction	 HANDS FREE MODE
Relaxing massage	Edematous cellulite	 BIPOLAR
<b>CLEAVAGE</b>	Fibrous cellulite	 COLOR DISPLAY
Face-lift	Sclerotic cellulite	 TOUCH SCREEN DISPLAY
Scars	Massage	 MAINS POWER SUPPLY
Firming	<b>BACK</b>	
<b>ARM-HAND</b>	Firming	
Wrinkles	Scars	
Firming	Massage	
Rejuvenation	Fat reduction	
Scars	<b>THIGHS AND LEGS</b>	
Stretch marks	Stretch marks	
Fat reduction	Firming	
<b>ABDOMEN-HIPS</b>	Scars	
Fat reduction	Lymphatic drainage	
Firming	Fat reduction	
Scars	Edematous cellulite	
Stretch marks	Fibrous cellulite	
	Sclerotic cellulite	
	Massage	

The **Pre-set/Aesthetic Area** includes 13 default programs intended for specific pathologies and connected with 5 different treatment categories.

TREATMENT CATEGORIES	PROGRAMS	TREATMENT CATEGORIES	PROGRAMS
<b>Lymphatic drainage</b>	1	<b>Lassità cutanea</b>	
<b>Localized fat</b>		• Minor	1
• Fold < 4 cm	1	• Moderate	1
• Fold 4-8 cm	1	• Serious	1
• Fold > 8 cm	1	<b>Face-cleavage rejuvenation</b>	
<b>Cellulite treatment</b>		• Minor	1
• Minor cellulite	1	• Moderate	1
• Moderate cellulite	1	• Serious	1
• Serious cellulite	1		

