

MAGNETOVET 200 with THERAPEUTIC BED

MAGNETOTHERAPY

MAGNETOVET 200 IS A HIGH INTENSITY PULSED ELECTROMAGNETIC FIELD MAGNETOTHERAPY EQUIPMENT. THE IMPULSES PRODUCED HAVE A THERAPEUTIC AND REGENERATIVE EFFECT ON THE MAIN OSTEOARTICULAR PATHOLOGIES THAT AFFECT PETS. FOR ITS SIMPLICITY OF USE MAGNETOVET 200 IS PARTICULARLY SUITABLE ALSO FOR HOME TREATMENT OF CHRONIC PAIN.

SPECIFIC PROGRAMS FOR DOGS, CATS AND HORSES

TECHNICAL FEATURES

Display

Backlit, 2.6"

Channels

2

Emitted intensity

200 average Gauss per channel with flex solenoids

Frequency

Adjustable in programming from 5 to 200Hz

Duty Cycle

From 5% to 50% in steps of 5

Power supply

- Mains
- Battery pack (optional)

Size


100x160x35 mm

Weight

440 gr

STANDARD EQUIPMENT

- 1 Carrying case
- 1 MAGNETOVET 200 device
- 1 Therapeutic bed 700x440 mm
- 1 Flexible solenoids 300x100 mm
- 1 Soft fleece cover
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTIONS

- Cod. G5620 Equipment without therapeutic bed
- Cod. G6039 Memory-foam therapeutic bed 4 coils, 700x440mm for cats and small dogs. For large dogs, buy 2 single beds and put them together with the velcro.



LAST 10



RENTAL OPTION



MULTI USER



1+1 MODE



PROGRAMMABLE



MAINS POWER SUPPLY

MAGNETOVET 200 with THERAPEUTIC BED

MAGNETOTHERAPY

Rev. En_03_2021

PROGRAM LIST

DOG / CAT

Fractures

Pseudoarthrosis

Wounds / Ulcers

Hip osteoarthritis

Knee osteoarthritis

Elbow osteoarthritis

Shoulder osteoarthritis

Spondylarthrosis

Luxation

Cicatrization

Edema

Peripheral paralyses

Inflammation pain

HORSE

Osteoarthritis

Wounds / Ulcers

Discopathy

Edema