

# MEDISOUND 3000

## ULTRASOUND THERAPY | 49 PROGRAMS

TOP PRODUCT SUITABLE FOR PHYSIOTHERAPY STUDIES AND HOSPITALS. ITS COMPACTNESS ALLOWS THE TRANSPORT IN ALL THE OPERATIONAL SITUATIONS.

31 BEAUTY | 18 HEALTH

### TECHNICAL FEATURES

#### Display

Visible area size 2.6"

#### Output

1

#### Hand

1 hand Ø 42 mm - IP68

#### Power

3 W/cm<sup>2</sup> ± 20%

#### Contact management

Yes

#### Emission Frequency

3 MHz and 1 MHz

#### Emission

Continuous and pulsed

#### Duty Cycle

10%-100%

#### Power supply

- Mains
- Rechargeable batteries (not included)

#### Size



170x220x60 mm

#### Weight

820 gr

### EQUIPMENT

- 1 Bag
- 1 MEDISOUND 3000 ultrasound therapy unit
- 1 Head
- 1 Gel flask
- 1 Power supply unit
- 1 Operating manual



CE  
0476



PROGRAMMABLE



MAINS SUPPLY



RECHARGEABLE BATTERIES

**GLOBUS**  
ITALIAN EXCELLENCE

# MEDISOUND 3000

## ULTRASOUND THERAPY | 49 PROGRAMS

### PROGRAMS LIST

#### HEALTH

Fracture

Wrist fracture

Tibia fracture

Carpal tunnel

Epicondylitis

Back pain

Calcific tendinitis

Muscle pain

Cervical pain

Trapezius pain

Spinal stenosis

Osteoarthritis

Spine osteoarthritis

Shoulder osteoarthritis

Ankle osteoarthritis

Knee osteoarthritis

Ankle osteoarthritis

Pseudoarthritis

#### BEAUTY

Oedematose cellulitis

Fibrous cellulitis

Lypolisis

Hypotonie

Vergetures

Skin orange

Skin smoothing

Lymphatic drainage

Massage

Skin restructure

Face lifting

Treatment photoaging

Abdomen loss of tone

Adipose abdomen

Culotte de cheval

Drainage thigh waist

Massage thigh waist

Adiposity buttocks waist

Relaxation

lymphatic drainage arms

Stretch mark breast

Breast loss of tone

Wrinkle treatment

Scare treatment

Couperose

Teleangectasie

Pimpled body

Pimpled face

Acne face

Cheloide treatment

Application active