

GENESY SII

ELECTROTHERAPY 2 CHANNELS | 60 PROGRAMS

PHYSIOTHERAPY TWO CHANNELS PRODUCT

16 TENS AND PAIN | 12 SPORT | 11 REHAB | 9 PREVENTION | 8 FITNESS, BEAUTY, WELLNESS
3 INCONTINENCE | 1 IONOPHORESIS

TECHNICAL FEATURES

Display

Visible area size 2,5"

Channels

2 independent (4 electrodes)

Frequency

1-150 Hz

Pulse width

50-400 μ s

Power

0-100 mA per channel

Power supply

Rechargeable battery

Size

80x130x25 mm

Weight

220 gr

EQUIPMENT

· 1 Bag

· 1 GENESY SII stimulator


· 2 Cables to connect the electrodes

· 4 Self-adhesive square electrodes

· 4 Self-adhesive rectangular electrodes

· 1 Charger

· 1 Operating Manual



CE
0476



IONO IONO

TENS TENS

PRO PROGRAMMABLE

Sys SYNCRO STIM

RECHARGEABLE BATTERIES

GLOBUS
ITALIAN EXCELLENCE

GENESY SII

ELECTROTHERAPY 2 CHANNELS | 60 PROGRAMS

PROGRAMS LIST

TENS AND PAIN

Conventional antalgic Tens

Endorphinic Tens

Rotator cuff tendinitis

Muscle pain

Knee osteoarthritis

Post-surgery pain

Muscle injuries

Cervical pain

Trapezius pain

Sciatica

Lombago

Epicondylitis

Scapulohumeral periarthritis

Carpal tunnel

Osteoarthritis

Menstrual pain

SPORT

Demo

Motor point pen

Warm-up lower limbs

Pre-competition warm-up lower limbs

Maximum strength lower limbs

Endurance strength lower limbs

Explosive strength lower limbs

Reactivity lower limbs

Capillarization lower limbs

Active recovery lower limbs

Muscle cool down lower limbs

Decontracting lower limbs

REHAB

Vastus medialis reinforcement

Swollen ankles

Hemiplegia upper limbs

Hemiplegia lower limbs

Recovery after ACL surgery

Functional recovery lower limbs

Quadriceps atrophy with prosthesis

Shoulder subluxation prevention

Muscular spasms

Multiple sclerosis - leg flexor extensor

Multiple sclerosis - ankle flexor extensor

PREVENTION

Ankle prevention

Knee prevention

Quadriceps muscle prevention

Reathletization lower limbs

Shoulder reathletization

Shoulder prevention

Elbow prevention

Reathletization upper limbs

Ankle muscle tone recovery

FITNESS, BEAUTY, WELLNESS

Leg anticellulite lipolysis

Abdominal area post-pregnancy lipolysis

Leg drainage

Abdominal area post-pregnancy drainage

Leg firming

Leg toning

Abdominal area post-pregnancy toning

Leg sculpting

INCONTINENCE

Mixed incontinence

Stress incontinence

URGE incontinence

IONOPHORESIS

Ionophoresis