

ELITE SII

2 CHANNELS ELECTROTHERAPY | 100 PROGRAMS

TWO CHANNELS PRODUCT FOR EVERY TREATMENT NEEDS.


44 FITNESS, BEAUTY, WELLNESS | 30 SPORT | 16 TENS | 6 PREVENTION | 3 REHAB | 1 INCONTINENCE

TECHNICAL FEATURES

Display
Visible area size 2,5"
Channel
2 independent (4 electrodes)
Frequency
1-150 Hz
Pulse amplitude
50-400 μ s
Power
0-100 mA for channel
Power supply
Rechargeable batteries
Size
80x130x25 mm
Weight
220 gr

EQUIPMENT

- 1 Bag
- 1 Stimulator ELITE SII
- 2 Cables for electrode connection
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 1 Charger
- 1 Operating manual



CE
0476



- TENS** TENS
- PRO** PROGRAMMABLE
- SyS** SYNCRO STIM
- RECHARGEABLE BATTERIES**

GLOBUS
ITALIAN EXCELLENCE

PROGRAMS LIST

SPORT

Demo
Motor point pen
Warm-up lower limbs
Pre-compet. warm-up lower limbs
Maximum strength lower limbs
Endurance strength lower limbs
Explosive strength lower limbs
Reactivity lower limbs
Capillarization lower limbs
Active recovery lower limbs
Muscle cool down lower limbs
Endurance lower limbs
Decontracting lower limbs
Abdominal maximum strength
Abdominal endurance strength
Abdominal active recovery
Abdominal cool down
Abdominal capillarization
Abdominal decontracting
Warm-up upper limbs
Pre-compet. warm-up upper limbs
Maximum strength upper limbs
Endurance strength upper limbs
Explosive strength upper limbs
Reactivity upper limbs
Capillarization upper limbs
Active Recovery upper limbs
Cool Down upper limbs
Endurance upper limbs
Decontracting upper limbs

TENS AND PAIN

Conventional antalgic tens
Endorphinic tens
Rotator cuff tendinitis
Muscle pain
Knee osteoarthritis
Post-surgical pain
Muscle injuries
Cervical pain
Trapezius pain
Sciatica
Lumbago
Epicondylitis
Scapulohumeral periarthritis
Carpal tunnel
Osteoarthritis
Menstrual pain

FITNESS, BEAUTY, WELLNESS

Leg anticellulite lipolysis
Hip-gluteus anticellulite lipolysis
Abdominal area anticellulite lipolysis
Upper limbs anticellulite lipolysis
Leg drainage
Hip-gluteus drainage
Abdominal area drainage
Drainage upper limbs
Leg firming
Hip-gluteus firming
Abdominal area firming
Firming upper limbs
Leg toning
Dorsal toning
Hip-gluteus toning
Abdominal area toning
Toning upper limbs
Leg sculpting
Hip-gluteus sculpting
Abdominal area sculpting
Leg vascularization
Vascularization upper limbs
Leg lipostress
Gluteus lipostress
Hip lipostress
Abdominal area lipostress
Lipostress upper limbs
Leg mass building
Mass building upper limbs
Abs area post- pregnancy drainage
Abs area post-pregnancy lipolysis
Abs area post- pregnancy toning
Man abdominal area definition
Man pectoral definition
Woman leg mass building
Woman upper limbs mass building
Swollen upper limbs
Swollen legs
Leg toning massage
Face beauty 1
Face beauty 2
Breast toning 1
Breast toning 2

PREVENTION

Ankle prevention
Knee prevention
Quadriceps muscle prevention
Shoulder reathletization
Shoulder prevention
Elbow prevention

REHAB

Quadriceps atrophy with prosthesis
Shoulder subluxation prevention
Ankle flexor-extensor

INCONTINENCE

Mixed Incontinence
