

# TECAR BEAUTY 7000 MED

## AESTHETIC TECARTHERAPY

THE BEST DEVICE WITH MEDICAL CERTIFICATION IN THE GLOBUS TECAR BEAUTY LINE. THE DEVICE IS EQUIPPED WITH A WIDE RANGE OF PROGRAMS INCLUDING BOTH AESTHETIC TREATMENTS AND PROTOCOLS FOR DIFFERENT MUSCULOSKELETAL PATHOLOGIES. THE TOUCHSCREEN NAVIGATION MENU ON A COLOR SCREEN INCLUDES AN ATLAS WITH SEVERAL PICTURES AND A PRACTICAL DESCRIPTION OF HOW TO PERFORM THE TREATMENTS. MOREOVER, TECAR BEAUTY 7000 MED ALLOWS EXPERT USERS TO CREATE AND STORE NEW PROGRAMS IN JUST A FEW STEPS.

### TECHNICAL SPECIFICATIONS

#### Display

Size 7" color touch screen

#### Nominal Power (NP)

350W ±10%

#### Emission frequencies

300kHz, 400kHz, 448kHz, 470kHz, 700kHz,  
1000kHz, 1200kHz, ± 10%

#### Treatment methods

Capacitive and Resistive

#### Capacitive Electrodes

3 Ø (30, 50, 70 mm)

#### Resistive Electrodes

3 Ø (30, 50, 70 mm)

#### Mains power supply

100-240V AC, 50-60Hz, 1.8A

#### Size

515x340x170 mm

#### Weight

11.5 Kg

### STANDARD EQUIPMENT

· Cod. G6244

· 1 TECAR BEAUTY BEAUTY 7000 MED device

· 1 Counter electrode (flexible plate)

· 1 Monopolar Resistive handpiece with 3 electrodes  
Ø (30, 50, 70 mm)

· 1 Monopolar Capacitive handpiece with 3 electrodes  
Ø (30, 50, 70 mm)

· 1 Bipolar handpiece Ø (30 mm)

· 1 Hands-free treatment kit

· 1 1000 ml conductive cream bottle

· 1 Power supply cable

· 1 Operating manual



# TECAR BEAUTY 7000 MED

## AESTHETIC TECARTHERAPY

### MEDICAL TREATMENT LIST

Rev. En\_03\_2024

Hallux valgus	Plantar fasciitis	Disc protrusion
Arthritis	Phlebitis	Acute osteitis pubis
Arthrosis	Fracture acute phase	Chronic osteitis pubis
Cervical arthrosis	Fracture sub-acute phase	Rhizartrrosis
Dorsal or lumbar arthrosis	Gonarthrosis	Subacromial impingement syndrome
Bursitis	1st degree muscular injury	Piriformis syndrome
Adhesive capsulitis	2nd degree muscular injury	Iliotibial band syndrome
Acute cervical pain	3rd degree muscular injury	Ulnar styloiditis
Chronic cervical pain	1st degree ligament injury	Acute Achilles tendinitis
Cervico-brachialgy	2nd degree ligament injury	Chronic Achilles tendinitis
Baker's cyst	3rd degree ligament injury	Acute tendinitis
Whiplash	Chronic lymphedema	Chronic tendinitis
Nerve compression	Post-traumatic lymphedema	Pes aserinus tendinitis
Patello-femoral chondropathy (cartilage lesion)	Acute Low Back Pain	De Quervain's tendinitis/Radial styloiditis
Muscle contracture	Chronic Lumbago	Patellar tendinitis
Coxarthrosis	Lumbosciatica	Tenosynovitis
Muscle cool down	Metatarsalgia	Sprain injuries with edema
Epicondylitis	Osteochondrosis/Osteochondritis	Sprain injuries without edema
Epitrochleitis	Periarthritis	Carpal tunnel
Disc herniation	Periostitis	
	Post-ACL surgery	















• The **Pre-set Area** includes 26 default programs intended for specific pathologies and connected with 5 different treatment categories.

TREATMENT CATEGORIES	SPECIFIC PATHOLOGIES	TREATMENT CATEGORIES	SPECIFIC PATHOLOGIES
<b>Lymphatic and vascular pathologies</b>		<b>Bone and articular pathologies</b>	
• Chronic lymphedema	1	• Arthritis	1
• Post-traumatic lymphedema	1	• Arthrosis	1
<b>Muscular and soft-tissue pathologies</b>		• Disc herniation	1
• Scars	1	• Fracture acute phase	1
• Muscle contracture	1	• Fracture sub-acute phase	1
• Muscle cool down	1	• Disc protrusion	1
• Fibrosis	1	• Sprain injuries	1
• 1st degree muscular injury	1	<b>Tendon and ligament injuries</b>	
• 2nd degree muscular injury	1	• Ligament injury 1st degree	1
• 3rd degree muscular injury	1	• Ligament injury 2nd degree	1
<b>Nervous pathologies</b>		• Ligament injury 3rd degree	1
• Cervico-brachialgy	1	• Acute tendinitis	1
• Nerve compression	1	• Chronic tendinitis	1
• Lumbosciatica	1		
• Piriformis syndrome	1		
• Carpal tunnel	1		

# TECAR BEAUTY 7000 MED

## AESTHETIC TECARTHERAPY

### ESTHETIC ATLAS PROGRAM LIST

<p><b>FACE</b></p> <p>Wrinkles</p> <hr/> <p>Eye contour</p> <hr/> <p>Acne scars</p> <hr/> <p>Lifting</p> <hr/> <p>Rosacea</p> <hr/> <p>Relaxing massage</p> <p><b>CLEAVAGE</b></p> <p>Face-lift</p> <hr/> <p>Scars</p> <hr/> <p>Firming</p> <p><b>ARM-HAND</b></p> <p>Wrinkles</p> <hr/> <p>Firming</p> <hr/> <p>Rejuvenation</p> <hr/> <p>Scars</p> <hr/> <p>Stretch marks</p> <hr/> <p>Fat reduction</p> <p><b>ABDOMEN-HIPS</b></p> <p>Fat reduction</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Stretch marks</p>	<p><b>GLUTEUS</b></p> <p>Stretch marks</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Lymphatic drainage</p> <hr/> <p>Fat reduction</p> <hr/> <p>Edematous cellulite</p> <hr/> <p>Fibrous cellulite</p> <hr/> <p>Sclerotic cellulite</p> <hr/> <p>Massage</p> <p><b>BACK</b></p> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Massage</p> <hr/> <p>Fat reduction</p> <p><b>THIGHS AND LEGS</b></p> <p>Stretch marks</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Lymphatic drainage</p> <hr/> <p>Fat reduction</p> <hr/> <p>Edematous cellulite</p> <hr/> <p>Fibrous cellulite</p> <hr/> <p>Sclerotic cellulite</p> <hr/> <p>Massage</p>	<p> TREATMENT ATLAS</p> <p> PULSED TREATMENTS</p> <p> ATHERMAL TREATMENTS</p> <p> RESISTIVE</p> <p> CAPACITIVE</p> <p> MULTI-FREQUENCY</p> <p> AUTOMATIC SCANNING OF FREQUENCIES</p> <p> MYOFASCIALS</p> <p> TRANSDERMAL DELIVERY</p> <p> HANDS FREE MODE</p> <p> BIPOLAR</p> <p> COLOR DISPLAY</p> <p> TOUCH SCREEN DISPLAY</p> <p> MAINS POWER SUPPLY</p>
---	--	--

The **Pre-set/Aesthetic Area** includes 13 default programs intended for specific pathologies and connected with 5 different treatment categories.

TREATMENT CATEGORIES	PROGRAMS	TREATMENT CATEGORIES	PROGRAMS
<b>Lymphatic drainage</b>	1	<b>Lassità cutanea</b>	
<b>Localized fat</b>		• Minor	1
• Fold < 4 cm	1	• Moderate	1
• Fold 4-8 cm	1	• Serious	1
• Fold > 8 cm	1	<b>Face-cleavage rejuvenation</b>	
<b>Cellulite treatment</b>		• Minor	1
• Minor cellulite	1	• Moderate	1
• Moderate cellulite	1	• Serious	1
• Serious cellulite	1		