

MAGNETOTHERAPY WELLNESS AT HOME

ARTHROSIS • JOINT PAIN
• OSTEOPOROSIS • BONE FRACTURES



GLOBUS
ITALIAN EXCELLENCE



SIMPLE AND EFFECTIVE MAGNETOTHERAPY IN THE COMFORT OF YOUR OWN HOME

HOW DOES IT WORK?

Our Magnum line device generate pulsed electromagnetic fields (PEMF) with low-frequency waves that can help bone and joint structures by strengthening their **microstructure when damaged** or **weak**.

Globus products are certified medical devices, **easy** to use and extremely **effective**.

IDEAL FOR

- Arthrosis and joint pain
- Osteoporosis
- Fractures and consolidation delays
- Bone oedema
- Inflammation

THE HEART OF GLOBUS PURE COPPER SOLENOIDS



More powerful
more effective
higher performance

An electromagnetic field is formed between the two pure copper components. The purer the copper, the more effectively and deeply the waves penetrate.

**This is Globus' secret:
the selection of the best components
and 100% Italian technology.**



CHOOSE YOUR GLOBUS HEART

The solenoid with the most powerful electromagnetic field on the market

Thanks to the double, fully conductive copper wire around the coils of our solenoids, Globus devices are able to emit a penetrative power with one of the highest levels on the market. The best way to treat illnesses.

FLEXIBLE SOLENOID



Very easy to position, it is ideal on joints (knees, ankles, shoulders, wrists, elbows) and small anatomical areas (feet and hands).

POCKET PRO SOLENOIDS



Soft and comfortable, adaptable to different anatomical sizes, they are ideal for treating the pelvis and back. Used in pairs, they generate a very large and penetrating electromagnetic field.

YOUR THERAPY AT YOUR FINGERTIPS

With **Globus Magnetotherapy**
it is even easier to feel well.
Effortlessly.

UP TO 70 PROGRAMS

Find your ideal device to treat
inflammation, arthrosis, joint pain,
fractures, bone oedema.

1+1 / 2+2 MODE

To launch 2 programmes at the same time
even with another person.

UP TO 4 OUTPUT CHANNELS

To connect different types of solenoids
and treat several body areas
simultaneously.

QUALITY ASSURANCE

MADE IN ITALY

ASSISTANCE FROM OUR SPECIALISED TECHNICIANS





MAGNUM 2 PRO DRIVE

Magnetotherapy is always with you.
Car, train, office: when and where you need it.

Magnum 2 Pro Drive is the most versatile and dynamic device in the Magnum range. Thanks to its built-in rechargeable battery, it can carry out therapy on the move. The device is fully equipped for extreme portability.


Its extensive range of preset programs makes it particularly effective against major musculoskeletal disorders.



IT'S ALL IN THE KIT

Your device always travels with you: a stylish and convenient case with everything you need for therapy on the go.

- Device
- Rechargeable battery
- "Soft 4" mat with 4 certified coils
- 2 "Flexible" solenoids
- USB cable for mobile charging
- Magnum test



CHOOSE YOUR IDEAL DEVICE

Buy or rent your magnetotherapy device



MAGNUM 3500 PRO

4

Output channels

Specifications

Intensity: 800 total average Gauss
Frequency: 5-200Hz

70

Programs

Supplied with

- Cod. G5434
- Device
- Case
- "Soft 4" mat
- 4 "Pocket Pro" solenoids
- 2 "Flexible" solenoids
- Magnum test
- Power supply unit
- Operating manual

2+2

2 + 2*
Mode

** Performs 2 different programs simultaneously, even on 2 people, with 2 channels each.*



MAGNUM 3000 PRO

2

Output channels

70

Programs

1+1

1+1* Mode

Specifications

Intensity: 400 total average Gauss
Frequency: 5-200Hz

Supplied with

- Device
- Bag
- Magnum test
- Power supply unit
- Operating manual

Purchase options

Code G5335

with 2 "Flexible" solenoids

Code G6034

with 2 "Pocket Pro" solenoids



Battery included



MAGNUM 2 PRO DRIVE

2 Output channels

41 Programs

1+1 1+1* Mode

Specifications

Intensity: 400 total average Gauss
Frequency: 5-200Hz

Supplied with

- Code G5429
- Device
- Bag
- Magnum test
- "Soft 4" mat
- 2 "Flexible" solenoids
- USB cable with 5-12 V connector
- Power supply unit
- Battery included
- Operating manual

MAGNUM XL PRO

2 Output channels

46 Programs

Specifications

Intensity: 400 total average Gauss
Frequency: 5-200Hz

Supplied with

- Device
- Bag
- Magnum test
- Power supply unit
- Operating manual

Purchase options

Code G3956
with 2 "Flexible" solenoids

Code G6033
with 2 "Pocket Pro" solenoids"

MAGNUM 2500

2 Output channels

52 Programs

1+1 1+1* Mode

Specifications

Intensity: 320 total average Gauss
Frequency: 5-200Hz

Supplied with

- Device
- Bag
- Magnum test
- Power supply unit
- Operating manual

Purchase options

Code G5438
with 2 "Flexible" solenoids

Code G6191
with 2 "Pocket Pro" solenoids

* Performs 2 different programs simultaneously, even on 2 people, with 2 channels each.



MAGNUM XL

2 Output channels

26 Programs

Specifications

Intensity: 280 total average Gauss
Frequency: 5-200Hz

Supplied with

Code G3216

- Device
- Bag
- "Flexible" solenoid
- Power supply unit
- Operating manual



MAGNUM L

1 Output channels

8 Programs

Specifications

Intensity: 120 total average Gauss
Frequency: 5-200Hz

Supplied with

Code G3947

- Device
- Bag
- "Flexible" solenoid
- Power supply unit
- Operating manual





TOTAL BODY MEMORY FOAM MATTRESS

Wellbeing as you sleep

If your therapy involves longer sessions or several body parts, GLOBUS has come up with a truly amazing **total body** solution.



SINGLE MAT

44x70 cm

To treat the spine and pelvis.

Code G5944
 "MAT 100 MEMORY FOAM"
 4 solenoids

MATTRESS

176x70 cm

For long-term treatments and treatments on several body areas.

Cod. G5947
 "TOTAL BODY 200 MEMORY FOAM" **Mattress**
 2 out of 4 active mats
 8 solenoids

Code G5948
 "TOTAL BODY 400 MEMORY FOAM" **Mattress**
 4 out of 4 active mats
 16 solenoids

Cod. G5949
 "TOTAL BODY 400 XP MEMORY FOAM" **Mattress**
 4 out of 4 active mats
 only compatible with
 MAGNUM 3500 PRO model



POSITIONING THE SOLENOIDS CORRECTLY MAKES THEM MORE EFFECTIVE



CERVICAL AREA



SHOULDER



UPPER BACK



ELBOW



LOWER BACK



HIP



KNEE



ANKLE



FOOT

ORIGINAL ACCESSORIES TO GUARANTEE WELLBEING



POCKET PRO SOLENOID

Code G6630

Very powerful and extremely comfortable solenoid, easy on joints and spine.

15,5x9,5 cm - 2 Coils



FLEXIBLE SOLENOID

Code G6627

Solenoid that simplifies application on joints (knees, ankles, shoulders, wrists, elbows) and small anatomical areas (feet and hands).

30x10 cm - 2 Coils



SOFT 4 SOLENOID

Code G6631

Accessory that makes therapy on the lower and upper back, neck and hips (pelvis) easier.

25x25 cm

4 Coils - 2 Channels



FLEXIBLE SOLENOID

Code G6629

Solenoid mainly suitable for spinal treatment.

40x15 cm - 4 Coils



FLEXIBLE SOLENOID COVER KIT

Code G5878

Practical sachets to keep the solenoids clean and ensure hygiene.



MAGNUM TEST

Code G6713

Checks the delivery of the magnetic field.

85x55 cm

YOU'LL FIND THE PROGRAM FOR YOU IN THE MAGNUM LINE

Discover all the illnesses in which magnetotherapy
can help you.



PROGRAMS	MAGNUM 3500 PRO 3000 PRO	MAGNUM 2 PRO DRIVE	MAGNUM XL PRO	MAGNUM 2500	MAGNUM XL	MAGNUM L
Algodystrophy / Sudeck's disease	✓	✓		✓		
Anti-inflammatory	✓	✓		✓		
Arthritis	✓	✓	✓	✓		✓
Arthrosis – Osteoarthritis	✓	✓	✓	✓	✓	✓
Ankle sprain	✓	✓	✓	✓		
Joint pain	✓	✓	✓	✓		
Bone oedem	✓	✓		✓		
Oedema-bruise	✓	✓	✓	✓		
Epicondylitis / Epitrochleitis	✓	✓	✓	✓		✓
Fracture with plaster	✓	✓	*	✓	*	*
Fractures	✓	✓	✓	✓	✓	✓
Spinal fusion	✓	✓				
Cartilage injury	✓	✓	✓	✓		
Backache	✓	✓	✓	✓	✓	✓
Osteonecrosis	✓	✓		✓		
Osteoporosis	✓	✓	✓	✓	✓	✓
Osteotomy	✓	✓				
Post-ACL surgery	✓	✓	✓	✓		
Hip-knee prostheses	✓	✓	✓	✓		
Pseudoarthrosis	✓	✓				
Delayed fracture consolidation	✓	✓	✓	✓	✓	
Shoulder tendonitis / Rotator cuff	✓	✓	✓	✓	✓	
Tendinitis	✓	✓	✓			
Carpal tunnel	✓	✓	✓	✓	✓	
Ulcers	✓	✓				
Multifrequency program	✓			1		
Preset frequency programs	25			15		
Night programs	3					
Cosmetic programs			14		10	

*** Use the “Fracture” program**

Some programs are arranged by body part. See the technical data sheets for more information.

GLOBUS

ITALIAN EXCELLENCE

[GLOBUSCORPORATION.COM](https://www.globuscorporation.com)



DOMINO S.R.L
Via Vittorio Veneto 52
31013 Codognè (TV) Italy

G6846 - REV. EN_04_2024

