

SHOCKWAVE 2000 - BeautyMed Line

RADIAL SHOCKWAVE DEVICE | 107 PROGRAMS

THE BEAUTYMED LINE SHOCKWAVE 2000 IS A GLOBUS RADIAL SHOCKWAVE DEVICE THAT EMITS HIGH-ENERGY ACOUSTIC PULSES USED IN THE DERMATOFUNCTIONAL AND MUSCULOSKELETAL FIELDS TO STIMULATE COLLAGEN PRODUCTION AND THE TISSUE REGENERATION PROCESS. IT IS A NON-INVASIVE TREATMENT METHOD THAT ACTS THROUGH MECHANICAL PRESSURE ON THE TISSUES. IN THE DERMATOFUNCTIONAL FIELD IT IMPROVES BLOOD CIRCULATION, HELPS TO REDUCE LOCALISED ADIPOSITY, TONES AND TIGHTENS THE SKIN. IN THE MUSCULOSKELETAL FIELD, IT IS PARTICULARLY EFFECTIVE IN TREATING TENDON AND MUSCLE PROBLEMS.

TECHNICAL FEATURES

Technology

Pneumatic

Pressure

1 - 5 Bar

Frequency

1 - 20Hz

Channels

2 (not usable simultaneously)

No. of programmable shocks

100-9900

Free Programs

User interface

- Color touchscreen tablet
- LCD screen on the handpiece

Handpiece included

Ergonomic, lightweight with integrated display and real time pressure indicator

Power supply

100-240V AC, 50-60Hz, 240VA

Dimensions LxPxH - Weight

430x320x150 mm - 16,6 Kg

STANDARD EQUIPMENT

- 1 ShockWAVE 2000 BeautyMed Line unit
- 1 Handpiece complete with projectile and guide tube
- 2 Handpiece holders
- 1 D36 head - Ø 36 mm
- 1 D20 head - Ø 20 mm
- 1 Tablet with software
- 1 Tablet holder
- 1 Spanner for handpiece maintenance
- 1 Cleaning brush for guide tube
- 1 Gel bottle 260 ml
- 1 Power cable
- 1 Operation guide
- 1 User manual

TRANSMITTERS HEADS SUPPLIED



D20 - Ø 20 mm



D36 - Ø 36 mm



0344



COLOR
DISPLAY



PROGRAMMABLE



TOUCH SCREEN
DISPLAY



RUN TIME



AUTO START



MAINS POWER
SUPPLY



GLOBUS
ITALIAN EXCELLENCE

SHOCKWAVE 2000 - BeautyMed Line

RADIAL SHOCKWAVE DEVICE | 107 PROGRAMS

41 DERMATOFUNCTIONAL PROGRAMS

FACE / DÉCOLLETÉ

Skin Rejuvenation

Acne scars

Elasticizing massage

Skin Tightening

Scar Treatment

UPPER LIMB

Elasticizing massage

Skin Tightening

Scar Treatment

Skin Rejuvenation

Fat Reduction

Lymphatic Drainage

HIPS

Lymphatic Drainage

Fat Reduction

Skin Tightening

Scar Treatment

Skin Rejuvenation

Elasticizing massage

ABDOMEN

Lymphatic Drainage

Fat Reduction

Skin Tightening

Scar Treatment

Skin Rejuvenation

Elasticizing massage

GLUTEUS

Lymphatic Drainage

Fat Reduction

Skin Tightening

Scar Treatment

Skin Rejuvenation

Elasticizing massage

Edematous cellulite

Fibrous cellulite

Sclerotic cellulite

THIGHS - LEGS

Fat Reduction

Skin Tightening

Scar Treatment

Skin Rejuvenation

Elasticizing massage

Edematous cellulite

Fibrous cellulite

Sclerotic cellulite

Lymphatic Drainage

SHOCKWAVE 2000 - BeautyMed Line

RADIAL SHOCKWAVE DEVICE | 107 PROGRAMS

66 MUSCULOSKELETAL PROGRAMS

NECK

Temporomandibular arthritis

Trigeminal neuralgia

Upper cervical trigger points

Lower cervical trigger points

SHOULDER

Frozen Shoulder

Cervicobrachial neuralgia

Dorsal scapular nerve entrapment

Hemiplegic shoulder pain

Shoulder myofascial pain syndrome

Nonunion of shoulder fracture

Rotator cuff calcification

Rotator cuff tendonitis (Acute)

Rotator cuff tendonitis (Chronic)

Subacromial impingement

Subacromial bursitis

Tenosynovitis of long head of biceps

ARM

Lateral epicondyle humerus fracture

Medial epicondyle humerus fracture

Lateral epicondylitis of the humerus

Medial epicondylitis of the humerus

FOREARM

Radial nonunion

Carpal tunnel syndrome

Nonunion of the wrist fracture

Cubital tunnel syndrome

HAND

Finger arthritis

De Quervain's Syndrome

Dupuytren disease

Trigger finger

BACK

Low Back Pain

Dorsal trigger points

Lumbar trigger points

HIP

Greater tuberosity fracture

Pressure ulcer

Bone marrow edema

Prostatitis

Femoral head aseptic necrosis

Gluteal myofascial pain syndrome

Bursitis greater trochanteric

Groin pain

Ischemic femoral necrosis

Pyriiformis syndrome

THIGH

Chronic proximal hamstring tendinopathy

Iliotibial band syndrome

Osteoporosis

Quadriceps tendonitis

Adductors trigger points

KNEE

Patellofemoral joint syndrome

Chronic synovitis of the knee

Delayed bone healing

Goosefoot bursitis

Meniscus tear

Myofascial pain syndrome of the knee

Osteoarthritis of the knee

Patellar tendinopathy

CALF

Calf myofascial pain syndrome

Tibial periostitis

Stress fracture of tibia

FOOT

Achilles tendinopathy

Diabetic foot ulcer

Heel spur

Metatarsal tendonitis

Morton's neuroma

Plantar fasciitis

Stress fracture

Anterior Talo-Fibular Ligament injury

Deltoid ligament injury

SHOCKWAVE 2000 - BeautyMed Line

RADIAL SHOCKWAVE DEVICE | 107 PROGRAMS

ACCESSORIES

THE TRANSMITTERS

Code **G7346**
LARGER TRANSMITTER D36
Diameter 36 mm



Code **G7347**
STANDARD TRANSMITTER D20
Diameter 20 mm



Code **G7348**
FOCUSED TRANSMITTER R15
Diameter 15 mm



Code **G7349**
SMALL TRANSMITTER D15
Diameter 15 mm



Code **G7350**
LONG TRANSMITTER F15
Diameter 15 mm



Code **G7351**
HIGH ENERGY D10
Diameter 10 mm



Code **G7352**
POINT TRANSMITTER A6
Diameter 6 mm



TO FACILITATE YOUR WORK

Code **G6510**
EASY GO CART



Code **G7378**
SHOCKWAVE TRANSPORT TROLLEY



Code **G7344**
ADDITIONAL HANDPIECE

