

SHOCKWAVE 2000 - Beauty Line

RADIAL SHOCKWAVE DEVICE | 17 PROGRAMS

THE SHOCKWAVE 2000 BEAUTY LINE IS A GLOBUS RADIAL SHOCKWAVE DEVICE THAT EMITS HIGH-ENERGY ACOUSTIC PULSES USED IN THE AESTHETIC FIELD TO STIMULATE COLLAGEN PRODUCTION AND TISSUE REGENERATION. IT IS A NON-INVASIVE TREATMENT METHOD THAT WORKS THROUGH MECHANICAL PRESSURE ON THE TISSUES, IMPROVING BLOOD CIRCULATION, HELPING REDUCE LOCALIZED FAT DEPOSITS, TONING, AND FIRMS THE SKIN.

TECHNICAL FEATURES

Technology

Pneumatic

Pressure

1 - 5 Bar

Frequency

1 - 20Hz

Channels

2 (not usable simultaneously)

No. of programmable shocks

100-9900

Free Programs

User interface

- Color touchscreen tablet
- LCD screen on the handpiece

Handpiece included

Ergonomic, lightweight with integrated display and real time pressure indicator

Power supply

100-240V AC, 50-60Hz, 240VA

Dimensions LxPxH - Weight

430x320x150 mm - 16.6 Kg

STANDARD EQUIPMENT

- 1 ShockWAVE 2000 Beauty Line unit
- 1 Handpiece complete with projectile and guide tube
- 2 Handpiece holders
- 1 D36 head - Ø 36 mm
- 1 D20 head - Ø 20 mm
- 1 Tablet with software
- 1 Tablet holder
- 1 Spanner for handpiece maintenance
- 1 Cleaning brush for guide tube
- 1 Gel bottle 1000 ml
- 1 Power cable
- 1 Operation guide
- 1 User manual

TRANSMITTERS HEADS SUPPLIED



D20 - Ø 20 mm



D36 - Ø 36 mm



COLOR
DISPLAY



PROGRAMMABLE



TOUCH SCREEN
DISPLAY



RUN TIME



AUTO START



MAINS POWER
SUPPLY



GLOBUS
ITALIAN EXCELLENCE

SHOCKWAVE 2000 - Beauty Line

RADIAL SHOCKWAVE DEVICE | 17 PROGRAMS

PROGRAM LIST

FACE / DÉCOLLETÉ

Deep oxygenation

Toning

Collagen activation

Rejuvenation

UPPER LIMB

Firming

Fat reduction

ABDOMEN

Firming

Slimming

HIPS

Decontracting

Firming

Slimming

GLUTEUS

Firming

Slimming

THIGHS / LEGS

Deep oxygenation

Firming

Heavy legs

Slimming

ACCESSORIES | THE TRANSMITTERS

Code **G7346**
LARGER TRANSMITTER D36
Diameter 36 mm



Code **G7347**
STANDARD TRANSMITTER D20
Diameter 20 mm



Code **G7348**
FOCUSED TRANSMITTER R15
Diameter 15 mm



Code **G7349**
SMALL TRANSMITTER D15
Diameter 15 mm



Code **G7350**
LONG TRANSMITTER F15
Diameter 15 mm



Code **G7351**
HIGH ENERGY D10
Diameter 10 mm



Code **G7352**
POINT TRANSMITTER A6
Diameter 6 mm



ACCESSORIES | TO FACILITATE YOUR WORK

Code **G6510**
EASY GO CART



Code **G7378**
SHOCKWAVES TRANSPORT TROLLEY



Code **G7344**
ADDITIONAL HANDPIECE

