

LASERVET 1.0

LASER THERAPY

LASERVET 1.0 IS THE ONLY LASER DEVICE IN CLASS 4 OF THE GLOBUSVET CATALOGUE. THE DEVICE CAN REACH 1 WATT OF POWER AND IT INCLUDES PROGRAMS FOR DOGS, CATS AND HORSES. MOREOVER, IT IS POSSIBLE TO CREATE AND SAVE NEW PROGRAMS.

SPECIFIC PROGRAMS FOR DOGS, CATS AND HORSES

TECHNICAL FEATURES

Display

Backlit, 2.6"

Laser source

Class 4 GaAlAs diode

Peak power

1 W $\pm 20\%$

Preset for wave length

808 nm

Issue

Pulsed and continuous

Target area dimensions

11 mm \varnothing , optics

Frequency

Up to 10000 Hz

Power supply

Mains and rechargeable battery

Size


100x160x35 mm

Weight

450 gr

STANDARD EQUIPMENT

- 1 Carrying case
- 1 LASERVET 1.0 device
- 1 Handpiece White 808 nm with 11 mm \varnothing , optics
- 1 Protective eyewear for the therapist
- 1 Kit with three small protective glasses for veterinary use, sizes S, M, L
- 1 Power supply unit
- 1 Operating manual



- LAST 10
- MULTI USER
- PROGRAMMABLE
- BATTERY PACK
- MAINS POWER SUPPLY

LASERVET 1.0

LASER THERAPY

MORE THAN 40 PROGRAMS FOR DOGS, CATS, HORSES.
TREATMENTS INCLUDE PROTOCOLS FOR:

Muscle-skeletal injuries such as arthrosis, tendinopathies, contusions, bursitis, sprains, muscle strains, etc.

Pathologies of the integumentary system such as skin or deep wounds

Pathologies of other systems such as otitis, sinusitis, stomatitis etc.

Rev. En_10_2025

