

LASERVET 15.0

LASER THERAPY

LASERVET 15.0 IS A PORTABLE DEVICE FOR HIGH PULSED POWER LASER THERAPY THAT ALLOWS TO PERFORM VETERINARY TREATMENTS IN QUICK TIMES. THE DISEASES TREATED ARE DIVIDED BY ANIMAL, WITH PARTICULAR ATTENTION TO EQUINES. THE ACTION OF LASER IS AIMED AT TISSUE REGENERATION, POST-TRAUMATIC AND ANTI-INFLAMMATORY TREATMENT AND REDUCTION OF BOTH CHRONIC AND ACUTE PAIN.

SPECIFIC PROGRAMS FOR DOGS, CATS, HORSES, SMALL MAMMALS, REPTILES AND BIRDS

TECHNICAL FEATURES

Display

Backlit, 2.6"

Laser source

Class 4 GaAlAs diode

Peak power

15 W $\pm 20\%$

Preset for wave length

808 nm

980 nm and 1064 nm (optional)

Issue

Pulsed and continuous

Target area dimensions

11 mm \varnothing and 30 mm \varnothing , optics

60 mm \varnothing , optics (optional)

Frequency

Up to 10000 Hz

Power supply

Mains and rechargeable battery


Size and weight

220x170x60 mm - 1.12 Kg



STANDARD EQUIPMENT

- 1 Carrying case
- 1 LASERVET 15.0 device
- 1 Handpiece HP 808 nm with 11 mm \varnothing and 30 mm \varnothing , optics
- 1 Protective eyewear for the therapist
- 1 Kit with three small protective glasses for veterinary use, sizes S, M, L
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTIONS

- Code **G6348** HP Handpiece 980 nm
- Cod. **G5996** Handpiece 1064 nm
- Cod. **G6089** Optics 60 mm \varnothing



LAST 10



MULTI USER



PROGRAMMABLE



BATTERY PACK



MAINS POWER SUPPLY

LASERVET 15.0

LASER THERAPY

Rev. En_10_2025

MORE THAN 130 PROGRAMS FOR DOGS, CATS, HORSES, SMALL MAMMALS, REPTILES AND BIRDS. TREATMENTS INCLUDE PROTOCOLS FOR:

Muscle-skeletal injuries such as arthrosis, tendinopathies, contusions, bursitis, sprains, muscle strains, fractures, dysplasia, neuropathies etc.

Pathologies of the integumentary system such as skin wounds, granulomas, sacculitis, perianal fistulas, pododermatitis, abrasions, etc.

Pathologies of other systems such as otitis, sinusitis, stomatitis etc.

