

MAGNETOVET 200 with THERAPEUTIC BED

MAGNETOTHERAPY

MAGNETOVET 200 IS A PORTABLE MAGNETOTHERAPY DEVICE CAPABLE OF GENERATING A HIGH-INTENSITY ELECTROMAGNETIC FIELD. THE PULSES IT PRODUCES HAVE A THERAPEUTIC AND REGENERATIVE EFFECT ON THE MAIN OSTEOARTICULAR CONDITIONS AFFECTING ANIMALS. THANKS TO ITS EASE OF USE, IN COMBINATION WITH THE THERAPEUTIC BED, MAGNETOVET 200 IS ALSO PARTICULARLY SUITABLE FOR HOME TREATMENTS OF CHRONIC PAIN.

SPECIFIC PROGRAMS FOR DOGS, CATS AND HORSES

TECHNICAL FEATURES

Display

Backlit, 2.6"

Channels

2

Emitted intensity

200 average Gauss per channel with flex solenoids

Frequency

Adjustable in programming from 5 to 200Hz

Duty Cycle

From 5% to 50% in steps of 5

Power supply

- Mains
- Battery pack

Size


100x160x35 mm

Weight

280 gr

STANDARD EQUIPMENT

- 1 Carrying case
- 1 MAGNETOVET 200 device
- 1 Therapeutic bed 700x440 mm
- 1 Flexible solenoids 300x100 mm
- 1 Soft fleece cover
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTIONS

- Cod. G5620 Equipment without therapeutic bed

Memory-foam therapeutic bed
4 coils, 700x440mm for cats and
small dogs.

- Cod. G6039 For large dogs, buy 2 single beds and
put them together with the velcro.



LAST 10



MULTI USER



PROGRAMMABLE



RENTAL
OPTION



1+1 MODE



MAINS POWER
SUPPLY



BATTERY PACK

MAGNETOVET 200 with THERAPEUTIC BED

MAGNETOTHERAPY

Rev. En_05_2024

PROGRAM LIST

DOG / CAT

Fractures

Pseudoarthrosis

Wounds / Ulcers

Hip osteoarthritis

Knee osteoarthritis

Elbow osteoarthritis

Shoulder osteoarthritis

Spondylarthrosis

Luxation

Cicatrization

Edema

Peripheral paralyses

Inflammation pain

HORSE

Osteoarthritis

Wounds / Ulcers

Discopathy

Edema